DAC-SpectroMed s.r.l.

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Cod 3004B100 100 ml
 Cod 3004B500 500 ml

PRINCIPLE

Direct bilirubin in the sample reacts with diazotized sulfanilic acid (DSA) forming a colored complex that can be measured by spectrophotometer. The intensity of coloration, measured at 540(±10) nm, is proportional to direct bilirubin.

CONTENTS AND COMPOSITION

Reagent AD
 Sulfanilic acid 30 mmol/l
 Hydrochloric acid 0.2 mol/l
 Reagent B
 Sodium nitrite 29 mmol/l

STORAGE AND STABILITY OF REAGENTS

Store reagents at 2-8°C and use until the expiry date shown on the label.

Symptoms of reagent unusability: presence of dredge, turbidity, blank absorption more than 0,03 at 540(±10) nm (1 cm cuvettes).

SAMPLES

Serum. Store protected from the light.
 Bilirubin in serum is stable within 3 months at -20°C, 4 days at 2-8°C, 1 day at 15-25°C.

REFERENCE VALUES

Adults³: Up to 0.2 mg/dl = 3.4 µmol/l.
 These values are approximate, it is recommended that each laboratory should establish its own reference values.

QUALITY CONTROL

It is recommended to use the Control Serum level I end II to verify the performance of the measurement procedure.
 Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT AND REAGENTS

Analyzer, spectrophotometer or photometer with filter 540(±10) nm
 Cuvettes with 1 cm light path. Thermostate for 37° C. Timer.
 Bilirubin Standard. Dropper for 50 µl, 100 µl and 1,5 ml.

PRECAUTION

The kit is only for in vitro use.
 All specimens must be considered potentially hazardous and handled as infectious.

REAGENT PREPARATION

Reagents are provided ready to use.
 Working Reagent preparation: 30 ml Reagent AD+1 ml Reagent B.
 Mix thoroughly. Stable for 10 days at 2-8°C.
 Sample Blank - Reagent **AD**.

NOTE

Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC cod 2051M5).

PROCEDURE

Method: end point
 Wavelength: 540(±10) nm
 Temperature: 37°C
 Cuvette: 1 cm light path
 Read against: Sample Blank

1. Pipette into labeled test tubes (Note):

	Sample Blank	Sample/ Standard
Reagent AD, ml	1,5	1,5
Sample/ Standard, µl	100	100
Reagent B, µl	-	50

NB: *Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.*

- Mix thoroughly and let the tubes stand for exactly 5 min at 37°C.
- Read the absorbance (As) of the Sample/Standard at 540(±10) nm against the Sample Blank.

CALCULATIONS

For the calculation of direct bilirubin concentrations in sample, calculate the difference between absorbencies of Sample (As) and Blank (Ab) for direct (DB) bilirubin.

Using these values, and calibration curve, plotted under total bilirubin standards, determine values of total and direct bilirubin.

SI UNITS: mg/dl direct bilirubin = µmol/l direct bilirubin : 17,1

METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0.05 mg/dl = 0.85 µmol/l

Linearity limit: 15 mg/dl = 257 µmol/l

At higher concentration dilute sample with distilled water 1:3 and repeat the measurement.

Repetition (within period):

Mean concentration	CV*	n*
0,59 mg/dl = 10,1 µmol/l	3%	20
6,74 mg/dl = 115,2 µmol/l	1%	20

Reproducibility (from period to period):

Mean concentration	CV*	n*
0,81 mg/dl = 13,8 µmol/l	7,1%	25
1,41 mg/dl = 24,1 µmol/l	6,3%	25

* CV - coefficient of variation n – number of determinations

Sensitivity: 48 mA . dl/mg = 2,81 mA l/µmol.

Interference:

Hemoglobin (10 g/l) does not interfere. Lipemia (triglycerides > 5 g/l) interfere. Other substances and drugs may interfere⁴.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Bilirubin is a waste product derived from the heme moiety of the hemoglobin released from senescent or damaged erythrocytes that are destroyed in the reticuloendothelial cells. After production, bilirubin is transported to the liver in association with albumin. Inside the hepatocytes bilirubin is conjugated with glucuronic acid and it is excreted into bile.

A number of inherited and acquired diseases affect production, uptake, metabolism, and excretion of bilirubin, resulting in hyperbilirubinemia^{3,5}.

Unconjugated hyperbilirubinemia is seen in newborns (physiological jaundice), in increased red cell destruction (hemolytic anemia, extensive hematoma), in ineffective erythropoiesis and in some rare genetic diseases (Gilbert's syndrome, Crigler-Najjar syndrome).

Conjugated hyperbilirubinemia is associated to a decreased excretion of bile due to liver diseases (hepatitis or cirrhosis) or to intrahepatic or extrahepatic cholestasis.

Jaundice is a clinical manifestation of hyperbilirubinemia, consisting of deposition of bile pigments in the skin, resulting in a yellowish staining of the skin and mucous membranes.

Clinical diagnosis must be made on integration of clinical and laboratory data.

BIBLIOGRAPHY

- Pearlman FC and Lee RTY. Detection and measurement of total bilirubin in serum, with use of surfactants as solubilizing agents. Clin Chem 1974; 20: 447-453.
- Zoppi F, Peracino A, Fenili D, Marcovina S and Ramella C. Metodo per la determinazione della bilirubina totale e coniugata. Uso di un tensioattivo cationico come agente solubilizzante. Giorn It Chim Cl 1976; 1:343-359.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.