
TRABECTOME®

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Trabectome Complications		
Total number of Patients	2454	100.0%
Trabeculectomy post Trabectome	153	6.2%
Shunt post Trabectome	63	2.6%
Diode Cyclo Photo Coagulation	7	0.3%
SLT post Trabectome	5	0.2%
Phaco post Trabectome	1	0.0%
ECP post Trabectome	9	0.4%
Trabectome post Trabectome (patient request)	22	0.9%
Tube Ligature Release on combined Tube and Trabectome	11	0.4%
Total Secondary Surgery	271	11.0%
Eyes out of study due to patient move - (no further follow up)	12	0.5%
Eyes out of study due to patients death (not related)	8	0.3%
Sustained Hypotony (IOP < 5 mmHg) one month post op	0	0.0%
Hypotony (IOP < 5 mmHg) one day post op	33	1.3%
IOP 10 > 10 mmHg from pre op IOP	143	5.8%
Aqueous Misdirection	1	0.0%
Infection	0	0.0%
Bleb Formation	0	0.0%
Wound Leaks	0	0.0%
Problematic Pain	0	0.0%
Choroidal Effusion	0	0.0%
Choroidal Hemorrhage	1	0.0%
Visual Acuity Decrease (> 2 lines)	0	0.0%
Intra Operative Blood Reflux		
Yes	1948	79.4%
No	80	3.3%
not reported	426	17.4%

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Mosaed S, Dustin L, Minckler DS. "Comparative Outcomes Between Newer and Older Surgeries for Glaucoma". Trans Am Ophthalmol Soc 2009;107:127-135

Minckler DS, Mosaed S, Dustin L, Francis BA and the Trabectome Study Group. "Trabectome (Trabeculotomy-internal approach): Additional experience and extended follow up." The American Ophthalmology Society, 2009

Liu J, Jung J, Francis BA. "Ab Interno Trabeculotomy: Trabectome Surgical Treatment for Open-angle Glaucoma." Expert Rev. Ophthalmol. 4(2), 119-128 (2009)

Minckler DS, Hill RA. "Use of Novel Devices for Control of Intraocular Pressure." Exp Eye Res. 2009 Apr;88(4):792-8. Epub 2008 Nov 30. Review

Patel SP, Sit AJ. "A Practical Model for Trabecular Meshwork Surgery." Arch Ophthalmol. 2009 Mar;127(3):311-3.

Godfrey DG, Fellman RL, Neelakantan A. "Canal Surgery in Adult Glaucomas" Curr Opin Ophthalmol. 2009 Mar;20(2):116-21. Review

Minckler DS, Mosaed S, Dustin L, Francis BA, Trabectome Study Group*. "Trabectome (Trabeculotomy-Internal Approach): Additional Experience and Extended Follow-Up". The American Ophthalmology Society, 2008;106:149-160

Nguyen QH. "Trabectome: a novel approach to angle surgery in the treatment of glaucoma." Int Ophthalmol Clin. 2008 Fa 11 ;48(4) :65-72. Review.

Francis BA, Minckler DS, Dustin L, Kawji S, Yeh J, Sit A, Mosaed S, Johnstone M, Trabectome Study Group*. "Combined Cataract Extraction and Trabectome Surgery (Trabeculotomy by Internal Approach) for Coexisting Cataract and Open-Angle Glaucoma". Journal of Cataract & Refractive Surgery, 2008; 34(7): 1093-1103

Filippopoulos T, Rhee DJ. "Novel Surgical Procedures in Glaucoma. Advances in Penetrating Glaucoma Surgery". Current Opinion in Ophthalmology, 2008; 19(2):149-54

Gunderson E, "Trabeculotomy Ab Interno, using the Trabectome: a promising treatment for patients with open-angle glaucoma." Insight. 2008 Jan-Mar;33(1):13-5

Mosaed S: "Ab Interno Trabeculotomy with the Trabectome Surgical Device". Techniques in Ophthalmology 5(2):63-66,2007

Francis BA, See RF, Rao N A, Minckler DS, Baerveldt G: "Ab-interno Trabeculectomy: Development of a Novel Device (Trabectome®) and Surgery for Open- Angle Glaucoma". Journal of Glaucoma, 2006; 15:68-73.

Minckler D, Baerveldt G, Ramirez MA, Mosaed S, Wilson R, Shaarawy T, Zack B, Dustin L, and Francis BA. "Clinical results with the Trabectome, a novel surgical device for the treatment of open-angle glaucoma." Trans Am Ophtlml Soc. 2006; 104:40-50.

Minckler DS, Baerveldt G, Ramirez-Alfaro M, Francis BA: "Clinical Results with the Trabectome" for Treatment of Open Angle Glaucoma". Ophthalmology, 2005; 112(6):962-967.

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Baerveldt G, Chuck R: "Minimally Invasive Glaucoma Surgical Instrument and Method". US Patent Number: US 6,979,328 B2 Issued: Dec. 27,2005.

Baerveldt G, Chuck R: "Minimally Invasive Glaucoma Surgical Instrument and Method". Canadian Patent Number: CA 2434173 Issued: Sept. 29,2009

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Trabectome® є зареєстрованою торговою маркою NeoMedix™ Corporation. © Корпорація NeoMedix™, 2010.