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PT MD 11-38623324-002:2002

## ACP - DAC

ACID PHOSPHATASE (TOTAL AND NON PROSTATIC)

KINETIC  $\alpha$ -NAPHTYL PHOSPHATE

For «in vitro» diagnostic use only  
Store at 2-8°C

Cod 2002A45 3x15 ml

### PRINCIPLE

Hillmann method: Phosphatase acid activity present in the sample is determined according to the modified method of Hillmann.

$\alpha$ -Nafthyl-phosphfate + H<sub>2</sub>O  $\xrightarrow{ACP}$   $\alpha$ -Naphtol + phosphate

$\alpha$ -Naphtol + Fast Red TR  $\longrightarrow$  Azo Dye

$\alpha$ -Naphtol reacts with a diazited compound forming a colour with a maximum of absorbance at 405(±10) nm.

Tartrate is used as specific of the prostatic fraction.

### CONTENTS AND COMPOSITION

Reagent 1	50 ml	pH 5,2
Sodium citrate buffer		50 mmol/l
Reagent 2	3 x 15 ml	
$\alpha$ -Nafthyl-phosphfate		10 mmol/l
Fast Red TR		6 mmol/l
Reagent 3	5 ml	
Sodium tartrate		2 mmol/l
Reagent 4	5 ml	
Acetic acid		0,5 mol/l

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

### Indications of deterioration:

Presence of particulate material, turbidity, absorbance of the blank  $\geq$  0,44 at 405(±10) nm (1 cm cuvette).

### SAMPLES

Serum<sup>1</sup>. Use only clear and unhemolyzed serum, seprated from the clotas soon as possible. Do not use plasma.

Acid phosphatase is very labile; stabilize by adding 50  $\mu$ l of Reagent 4 per 1 ml of the sample. Stability for 7 days at 2-8°C.

### REFERENCE VALUES<sup>4,5</sup>

	30°C	37°C
Acid phosphatase Total		
Men:	< 4,3 U/l	< 5,4 U/l
Women:	< 3,1 U/l	< 4,2 U/l
Acid phosphatase Prostatic	< 1,5 U/l	< 1,7 U/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 405(±10).  
Pipette for 10  $\mu$ l, 100  $\mu$ l, 1,0 ml.

### PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious.

### REAGENT PREPARATION

Working Reagent: Dissolve one tablet of Reagent 2 in 15 ml of Reagent 1. Cap vial and mix gently to dissolve contents.

Stability of Working Reagent – 2 daus at 2-8°C.

### PROCEDURE

Assay conditions

Method:	kinetic
Wavelength :	405 (±10) nm
Light path:	1 cm
Temperature :	30°C/37°C
Blank:	air or distillate water

1. Bring the Working Reagent and the instrument to reaction temperature.

2. Pipette into a cuvette:

Acid phosphatase	Total	Non Prostatic
Working Reagent	1,0 ml	1,0 ml
Reagent 3	-	10 $\mu$ l
Sample	100 $\mu$ l	100 $\mu$ l

NB: Volumes of reagents and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix thoroughly and let stand the tubes for 5 minute.

4. Read initial absorbance (A) of the sample, start the stopwatch and read absorbance at 1 minute intervals thereafter for 3 minutes.

5. Calculate the difference between absorbance and the average absorbance difference per minute ( $\Delta A/min$ ).

### CALCULATIONS

The acid phosphatase concentration in the sample (U/l) is calculated using the following general formula:

$$C \text{ Total} = \Delta A / \text{min Total} \times 750$$

$$C \text{ Prostatic} = (\Delta A / \text{min Total} - \Delta A / \text{min Non prostatic}) \times 750$$

750 – Factor 1 U/l = 16,67 nkat/l.

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,13 U/l.

Linearity limit: 150 U/l.

Repeatability (within run):

Mean Concentration	CV	n
23,67 U/l	0,95 %	20
2,56 U/l	2,9 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
23,6 U/l	0,92 %	20
2,6 U/l	2,76 %	20

CV - coefficient of variation n – number of determinations

Sensitivity: 1 U/l = 0,0034 A/min.

Interferences: Hemolysis interferes due the high concentration of acid phosphatase in red cells<sup>1</sup>. A list of drugs and other interfering substances with acid phosphatase determination has been reported by Young et. al<sup>2,3</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Acid phpsphatase is an enzyme present in almost all weaves of the organism, being particulary high in prostate, stomach, liver, muscle, spleen, erythrocytes and platelets. High levels of acid phosphatase are found in prostatic patologies as hypertrophy, prostatitis or carcinoma. In hematological disorders, bones or liver diseases as well as in Paget's or Gaucher's diseases.

Decreased serum acid phosphatase has no clinical significance<sup>1,4,5</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. Abbot L, et al. Acid phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co, St Louis. Toronto. Princeton 1984; 1079-1083.
2. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCC Press, 1995.
3. Young DS. Effects of disease on clinical laboratory tests, 4th ed. AACCC, 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd edition. AACCC 1999.
5. Tietz et al. Clinical Guide to Laboratory Tests, 3rd. AACCC 1999.
6. J.Henry Wilkinson. The principles and practice of Diagnostic Enzymology. Edward Arnold ed., 1981.