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PT MD 11-38623324-002:2002

ALP-DEA-DAC

ALKALINE PHOSPHATASE
 KINETIC TEST WITH DIETHANOLAMINE BUFFER

For «in vitro» diagnostic use only
 Store at 2-8°C

Cod 2006A100 100 ml

PRINCIPLE

In alkaline medium alkaline phosphatase (ALP) catalyzes the transfer of the phosphate group from 4-nitrophenylphosphate to diethanolamine (DEA), releasing 4-nitrophenol. The activity of ALP is determined from the rate of colour formation, measured at 405-410 nm^{1,2}.

CONTENTS AND COMPOSITION

Reagent A	1 x 80 ml	pH 9.8
Diethanolamine		1,0 mol/l
Magnesium chloride		0,6 mmol/l
Reagent B	1x 20 ml	
4-nitrophenylphosphate		10 mmol/l

Harmful! Avoid contact with skin and mucous membranes.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label. Indications of deterioration: presence of particulate material, turbidity, absorbance of working reagent ≥ 0.850 at 405 (± 10) nm (1 cm light path cuvettes).

SAMPLES

Serum free of Hemolysis.
 Alkaline Phosphatase in serum or plasma is stable for 7 days at 2-8°C.

REFERENCE VALUES

Children: 180-1200 U/l.
 Adult: 98-279 U/l.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 405-410 nm.
 Cuvettes with 1 cm light path.
 Pipettes for 20 μ l, 200 μ l and 800 μ l.

PRECAUTION

For in vitro diagnostics only.
 All patients' samples should be considered as potentially dangerous and handled as infectious. Follow safety measures for handling harmful and poisonous substances.

REAGENT PREPARATION

Reagents are ready to use.

PROCEDURE

Method: kinetic (increasing)
 Wavelength: 405-410 nm
 Temperature: 37°C
 Read against: distilled water or air

1. Bring the Reagents and the instrument to reaction temperature (37°C).

2. Pipette into a cuvette 1 cm light path:

Reagent A 800 μ l

Sample, Standard 20 μ l

3. Mix and incubate for 1 minute. Add in cuvette:

Reagent B 200 μ l

4. After 1 minute record the initial absorbance against distilled water, then measure absorbance for 2 minutes at 1 minute interval.

5. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ($\Delta A/\text{min}$).

*NB: reagents and sample volumes can be changed proportionally according to the working volume of the cuvette.

CALCULATIONS

The ALP concentration in the sample (U/l) is calculated using factor:

$$405 \text{ nm: Activity (U/l)} = \Delta A / \text{min}_{\text{Sam}} \times 3250$$

METROLOGICAL CHARACTERISTICS

Detection limit: 0,001 $\Delta A/\text{min} = 3,25$ U/l.

Linearity limit: 0,554 $\Delta A/\text{min} = 1800$ U/l.

Repeatability (within run):

Mean Concentration	CV	n
191 U/l	0,91 %	20
642 U/l	0,96 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
176 U/l	2,29 %	25
416 U/l	1,43 %	25

CV - coefficient of variation n - number of determinations

Interferences: Bilirubin up to 128,3 $\mu\text{mol/l}$ (0,075 g/l), lipids up to 10 g/l, glucose up to 55,5 mmol/l (10 g/l) and ascorbic acid up to 2,84 mmol/l (0,5 g/l) do not interfere. Other drugs and substances may interfere³.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Alkaline Phosphatase catalyzes the hydrolysis of organic phosphate monoesters at alkaline pH. The enzyme is present in practically all tissues of the body, especially at or in the cell membranes, and it occurs at particularly high concentrations in placenta, intestinal epithelium, kidney tubules, osteoblasts and liver. The form present in the sera of normal adults originates mainly in the liver and bone.

Elevated serum ALP is found in patients with bone disease associated with increased osteoblastic activity (Paget's disease, primary and secondary hyperparathyroidism, bone tumors, rickets, osteomalacia, bone fractures) and also in patients with hepatobiliary disease (obstructive jaundice, hepatitis, hepatotoxicity caused by drugs, liver cancer). Physiological changes, such as bone growth and pregnancy, may cause increases in ALP levels^{5,6}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Haussament T.U. et. Al. Clin. Chem. Acta 35. 271-273 (1977).