DAC-SpectroMed s.r.l.

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PT MD 11-38623324-002:2002

Amylo-DAC

α -AMYLASE FOR SERUM AND URINE
AMYLOCLASTIC METHOD BY KARAVEY

For «in vitro» diagnostic use only
Store at 2-8°C

Cod 2015A200 200 ml

PRINCIPLE

α -amylase catalyzes hydrolysis of starch and forms end products that do not form coloration with iodine.

Excess of starch reacts with iodine and form colored complex.

α -amylase activity is judge by decreasing of intensity of coloration, measured at 630-670 nm.

CONTENTS AND COMPOSITION

Reagent A	40 mg
Water-soluble starch	
Reagent B	7,56 g, при растворении:
Sodium phosphate	0,2 mol/l
Benzoic acid	0,07 mol/l
Reagent C	10 ml
Iodide	0,1 n

STORAGE AND STABILITY OF REAGENTS

Reagents are stable up to the indicated expiry date if stored at 2-8°C.

Substrate-buffer solution is stable 1 month at 18-22°C.

Working Reagent C is stable 24 hours at 18-22°C.

SAMPLES

Serum, plasma and urine.

α -amylase in serum, plasma and urine is stable 5 days at 2- 8°C.

Heparin is recommended as anticoagulant.

REFERENCE VALUES

Serum: 16 - 30 g/l-hour Urine: 28 - 160 g/l-hour

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer, photometer able to read at 630-670 nm. Thermostat at 37°C.

Pipette for 10 μ l, 0,50 ml and 4,0 ml.

PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Substrate-buffer solution:

Content of Reagent B bottle transfer quantitatively in the thermostable flask and add 80 ml of distilled water.

Dissolve under heating and bring to boil.

Content of Reagent A bottle mix with 3-5 ml of cool distilled water and insert in boiling Reagent B under constant mixing.

Boil during 1 minute, cool off, add distilled water up to 100 ml and mixed well. Store tightly closed.

Working Reagent C:

To dissolve Reagent C in distilled water in 10 times:

1 ml Reagent C + 9 ml of distilled water.

PROCEDURE

Assay conditions

Method:	end point
Wavelength :	630-670 nm
Light path:	1 cm
Temperature :	37°C
Blank:	distilled water

1. Pipette into labeled test tubes:

	Control	Test
Substrate-buffer solution, ml	0.50	0.50

2. Incubate 5 minutes at 37°C in thermostat.

Pipette into labeled test tubes:

	Control	Test
Sample, ml	-	0,01

4. Start the stopwatch in the moment of sample inserting in substrate-buffer solution. Incubate exactly 5 minutes at 37°C in thermostat.

5. Pipette into labeled test tubes:

	Control	Test
Working Reagent C, ml	0.5	0.5
Distilled water, ml	4.0	4.0

6. Measured absorbencies of the Test and Control at 630-670 nm against distilled water.

If cuvettes with another working volume are used recommended proportionally change volumes of components.

CALCULATIONS

α -amylase activity in the sample in 1g of starch hydrolyzed with 1 l of serum during 1 hour of incubation at 37(\pm 1)^oC calculated as a following:

$$K_s = \frac{A_c - A_t}{A_c} \times 240$$

K_s - α -amylase activity in the sample, g/l-hour;

A_c – absorbance of the Control;

A_t – absorbance of the Test;

240 = 0,2 x 1200

0,2 – amount of starch (mg), inserting in the Control and in the Test;

1200 = 100 x 12

100 – coefficient of recalculation from 0.2 mg to 1g of starch hydrolyzed with 1 l of serum

12 – coefficient of recalculation from 5 minutes to 1 hour.

$$1 \text{ g/l } \times \text{ h} = 1 \text{ mg/ml } \times \text{ h} = 0,278 \text{ mg/l } \times \text{ c};$$

$$1 \text{ mg/l } \times \text{ c} = 3,6 \text{ mg/ml } \times \text{ h} = 3,6 \text{ mg/ml } \times \text{ h}$$

METROLOGICAL CHARACTERISTICS

Detection limit: 10,8 g/l x h.

Linearity limit: 130 g/l x h. For higher values dilute sample

with distilled water and repeat measurement.

CV: < 7 %.

Interferences:

Hemolysis doesn't interfere.

Saliva and skin content a-amylase.

Avoid contacts with skin and saliva³.

DIAGNOSTIC CHARACTERISTICS

δ -amylase catalyzes hydrolysis of δ 1,4-linkage carbohydrates which consist of δ -D-glucose molecules. As a result, dextrans, maltose and glucose molecules are formed. δ -amylase is secreted by exocrine part of pancreas (P-type) and by salivary glands (S-type), but also is found in other tissues of organism.

Measurements of amylase are used primarily in the diagnosis and treatment of the diseases of the pancreas. Amylase is found primarily in the pancreas and salivary glands. When released in the digestive tract, the enzyme hydrolyzes starch.

Amylase determinations are useful in the diagnosis of diseases of the pancreas and parotids. Elevated serum levels are associated with acute pancreatitis and other pancreatic disorders as well as mumps and bacterial parotitis. Hyperamilasemia may be caused by renal insufficiency, acute pain in abdominal cavity, lungs and ovary tumors, salivary glands injury, macroamilasemia, ketoacidosis, bile-duct disease, brain traumas, chronic alcoholism and consumption of opiates^{4,5}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

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