DAC-SpectroMed s.r.l.

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
PRINCIPLE

δ -Amylase catalyses hydrolysis 2-chlor-4-nitrophenilmaltotriozide (CNP-G3) with forming of 2-chlor-4-nitrophenol (CNP). The intensity of CNP coloration, measured at 405(\pm 10) nm, is proportional to δ -amylase activity.

CONTENTS AND COMPOSITION

Reagent	pH 6,00
CNP- G3	4,55 mmol/l
Calcium acetate	5 mmol/l
Sodium chloride	51,5 mmol/l
Buffer	50 mmol/l
Preservatives	

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

Indications of deterioration:

Reagents: presence of particulate material, turbidity, absorbance of the blank over 0,50 at 405(\pm 10) nm (1 cm cuvette).

SAMPLES

Serum, urine.

α -amylase in serum, plasma and urine is stable 1 days at 2- 8°C. Heparin is recommended as anticoagulant.

REFERENCE VALUES

Serum: 15-100 U/l

Urine: \leq 400 U/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 405(\pm 10) nm. Cuvettes with 1 cm light path. Pipettes for 10 μ l and 600 μ l.

PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious.

REAGENT PREPARATION

The reagent is ready to use.

PROCEDURE

Assay conditions

Method:	kinetic (increasing)
Wavelength :	405 (\pm 10) nm
Light path:	1 cm
Temperature :	37°C
Blank:	distilled water

1. Bring Reagent and photometer to reaction temperature.
2. Pipette into a cuvette:

Reagent	600 μ l
Sample	10 μ l

NB: Volumes of reagent and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix thoroughly and place cuvette into photometer. Start timer.
4. After 1 minute record initial absorbance against the distilled water and record absorbance at 30 sec, intervals thereafter for 2 minutes.
5. Calculate the difference between consecutive absorbances and mean absorbance difference for 1 minute ($\Delta A/\text{min}$).

CALCULATIONS

δ -Amylase activity in sample is calculated by formula:

$$\frac{\Delta A/\text{min}_{\text{sam}}}{\Delta A/\text{min}_{\text{st}}} \times C_{\text{st}} = C_{\text{sam}}$$

Calculation using factor

$$405 \text{ nm: Activity (U/l)} = \Delta A/\text{min}_{\text{sam}} \times 4325$$

METROLOGICAL CHARACTERISTICS

Detection limit: 0,001 $\Delta A/\text{min}$ = 4,325 U/l.

Linearity limit: 0,694 $\Delta A/\text{min}$ = 3000 U/l.

Coefficient of variation: < 2 %.

Interference:

Bilirubin (60 mg/dl), lipid (10 g/l), glucose (20 g/l) and ascorbic acid (1 g/l) do not interfere. Other drugs and substances may interfere.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

δ -amylase catalyzes hydrolysis of δ 1,4-linkage carbohydrates which consist of δ -D-glucose molecules. As a result, dextrans, maltose and glucose molecules are formed. δ -amylase is secreted by exocrine part of pancreas (P-type) and by salivary glands (S-type), but also is found in other tissues of organism.

Appreciation of δ -amylase activity in serum and urine is widely used for diagnostics of pancreas disease, such as acute and chronic pancreatitis. Hyperamilasemia may be caused by renal insufficiency, acute pain in abdominal cavity, lungs and ovary tumors, salivary glands injury, macroamilasemia, ketoacidosis, bile-duct disease, brain traumas, chronic alcoholism and consumption of opiates.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

Chelating agents interfere with the reaction. Do not use citrate, oxalate or EDTA anti-coagulant.

The reagent contains calcium, which can cause the precipitation of the fibrinogen from plasma.

BIBLIOGRAPHY

1. Street, H.V. and Close, J.R., Clin. Chem. Acta 1:256 (1956).
2. Henry, R.J. and Chiamori, N. Clin. Chem. 6:434 (1960).
3. David, H., Clin. Chem. 28:1485 (1985).
4. McCroskey, R., Chang, T., David, H. and Winn,E., Clin. Chem. 28:1787, (1982).
5. Young DS. Effects of drugs on clinical laboratory tests 3rd ed. AACC Press, 1990.
6. Balsells D, Gella FJ, Gubern G, Canalias F. Reference values for α amylase in human serum and urine using 2-chloro-4-nitrophenyl- α -D-maltotriozide as substrate. Clin Chim Acta 1998; 274: 213-217.
7. Young DS. Effects of drugs on clinical laboratory tests 4th ed. AACC Press, 1995.
8. Tietz Textbook of Clinical Chemistry 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
9. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.