DAC-SpectroMed s.r.l.

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PT MD 11-38623324-002:2002

Amylo Pancrea-DAC.Lq

PANCREATIC α -AMYLASE IN SERUM AND URINE
 KINETIC TEST WITH EPS-G7 SUBSTRATE

For «in vitro» diagnostic use only
 Store at 2-8°C

Code 2021A75 75 ml

METHOD PRINCIPLE

Enzymatic photometric test where 4,6-ethylidene-(G7)-p-nitrophenyl-(G1)- α -D-maltoheptaoside (EPS-G7) substrate is decomposed into different fragments by α -Amylase. These fragments are further hydrolyzed by α -glucosidase producing glucose and p-nitrophenol. As salivary isoenzyme is selectively inhibited by two monoclonal antibodies combination before incubation, the increase of absorption expresses pancreatic amylase activity in the sample.

CONTENTS AND COMPOSITION

Code		2021A75
Reagent A	pH=7,15	3x20 ml
Good's buffer	0.1 mol/l	
Sodium chloride	62.5 mmol/l	
Magnesium chloride	12.5 mmol/l	
Glucose oxidase	≥ 2 KU/l	
Monoclonal antibodies	≥ 25 mg/l	
Reagent B	pH=7,15	1x15 ml
Good's buffer	0.1 mol/l	
EPS-G7	8.5 mmol/l	

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

SAMPLES

Serum, plasma with heparin or EDTA, urine.
 α -amylase in serum and urine is stable 7 days at 2- 8°C.

REFERENCE VALUES

	Women	Men
Serum/plasma	< 53 U/l	< 53 U/l
Urine	< 319 U/l	< 356 U/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.
 Each laboratory should establish its own internal Quality Control scheme and procedures.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 405 (400-420) nm.
 Pipettes for 10, 20, 250 and 1000 μ l.

PRECAUTION

For in vitro diagnostics only.
 Handle all patients' samples as potentially dangerous and treat as infectious.

REAGENT PREPARATION

The reagents are ready to use.

PROCEDURE

Method: kinetic
 Wavelength : 405 nm
 Temperature : 37°C
 Blank: reagent

1. Pipette into labeled tubes:

	Blank	Serum/plasma	Urine
Sample	-	20 μ l	10 μ l
Reagent A	1000 μ l	1000 μ l	1000 μ l

2. Mix the contents, incubate for 3 minutes, then add:

Reagent A	250 μ l	250 μ l	250 μ l
-----------	-------------	-------------	-------------

NB: Volumes of reagent and samples can be proportionally changed according to the cells working volumes of used analyzers.

3. Mix thoroughly and incubate for 2 minutes at 37°C.

4. Read initial absorbance against blank, then read at 1 minute interval within 3 minutes.

5. Calculate the difference between consecutive absorbances and mean absorbance difference per minute ($\Delta A/min$).

CALCULATIONS

α -Amylase activity in sample is calculated by formula:

$$\text{Activity (U/l)} = \Delta A/\text{min} \times 5670$$

METROLOGICAL CHARACTERISTICS

Detection limit: 5 U/l

Linearity limit: $\Delta A/min = 0.350$. If this value is exceeded, sample should be diluted 1:10 by 9 g/l NaCl solution. Result should be multiplied by 11.

Intra-assay reproducibility:

Mean concentration	CV*	n*
69.7 U/l	3.13	20
207 U/l	1.26	20
370 U/l	0.91	20

Inter-assay reproducibility:

Mean concentration	CV*	n*
68.3 U/l	2.17	20
204 U/l	0.79	20
371 U/l	0.85	20

CV - coefficient of variation; n - number of tests

Interference: bilirubin (up to 40 mg/dl), lipemia (triglycerides up to 2000 mg/dl), and ascorbic acid (up to 30 mg/dl) do not interfere. Haemoglobin may interfere at minimum concentrations. Other drugs and substances may interfere.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Measurements of amylase activity are primarily used for diagnostics and treatment of pancreas disease.

Amylase is found primarily in pancreas and salivary glands. When released in the digestive tract this enzyme hydrolyzes starch. Amylase determinations are useful in the diagnosis of diseases of pancreas and parotids.

Elevated serum levels are associated with acute pancreatitis and other pancreatic disorder as well as mumps and bacterial parotitis. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

Residual activity of salivary α -amylase is up to 3%. Elevated levels of salivary α -amylase may lead to elevated results of pancreatic α -amylase activity. Saliva and skin contain α -amylase, that's why reagents should not come into contact with skin and be pipette by mouth.

BIBLIOGRAPHY

- Lorentz K α -Amylase. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998, p.192-202.
- Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood AR, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company; 1999, p.689-98.
- Gerber M, Naujoks K, Lenz H, Wulff K. A monoclonal antibody that specifically inhibits human salivary alpha-amylase. Clin. Chem. 1987; 33: 1158-62.
- Kruse-Jarres JD, Kaiser C, Halfkenschied JC, Hohenwallner W, Stein W., Bohner J et al. Evaluation of a new alpha-amylase assay using 4,6-ethylidene-(G7)-1-4-nitrophenyl-(G1)- α -D-maltoheptaoside as substrate. J Clin Chem Biochem 1989;27:103-13.
- Tietz NW, Burlina A, Gerhardt W, Junge W, Maffertheimer P, tvaluation of a specific pancreatic isoamylase assay based on a double monoclonal-antibody technique. Clin Chem 1988; 34:2096-102.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001. P 16-17, 50-51.