**DAC-SpectroMed s.r.l.**  
MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64  
Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920  
Email: office@dacspectromed.com  
www.dacspectromed.com  
PT MD 11-38623324-002:2002

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
Cod 2023A1200 400 tests

### PRINCIPLE

Aspartate aminotransferase (AST) catalyzes the reversible transfer of the aspartate amino groups to  $\alpha$ -ketoglutarate by reactions, described below. The activity of AST is proportional to content of dinitrophenylhydrazosone, measured at 500-550 nm.

### CONTENTS AND COMPOSITION

Reagent A	<b>pH 7,4</b>
L-aspartate	0,1 mol/l
$\alpha$ -ketoglutaric acid	2 mmol/l
phosphate buffer	0,1 mol/l
Reagent B	
2,4-dinitrophenylhydrazosone	1 mmol/l
Reagent C	
sodium hydroxide	2,0 mol/l
Reagent S	2,5 ml
calibration solution of sodium pyruvate	2,0 mmol/l

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

### SAMPLES

Serum free of Hemolysis.

AST in serum is stable for 7 days at 2-8°C.

### REFERENCE VALUES

0,10-0,45 mmol/h x l = 28-125 nmol/s x l = 1,7 - 7,5 IU/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with filter 500-550 nm. Water bath 37°C. Stopwatch.

Pipettes for 25,0  $\mu$ l, 0,25 ml, and 2,5 ml.

### PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

### REAGENT PREPARATION

Reagent A, Reagent B and Reagent S are ready to use.

Preparation of 0,4 M NaOH: Dilute Reagent C in volumetric flask with distilled water, free from carbonates, up to volume specified on a label. Store in polyethylene bottle and use until appearance of

### PROCEDURE

Testing of parallel samples is obligatory!

Assay conditions

Wavelength: 500-550 nm  
Temperature: 18-25/37°C  
Cuvette: 1 cm light path  
Read against: against reagent  
Method: Raitman-Frenkel

1. Pipette in labeled tubes:

	Tested sample	Control sample
Reagent A	0,25 ml	0,25 ml
Distilled water	-	0,05 ml
Blood serum	0,05 ml	-

2. Mix and incubate for 60 minutes on water bath at 37°C.

3. Add in labeled tubes 0,25 ml of Reagent B.

4. Mix and incubate for 20 minutes at 18-25°C.

5. Add in labeled tubes 2,5 ml of 0,4 M NaOH.

6. Mix and incubate for 10 minutes at 18-25°C.

7. Measure the absorbance of tested sample relative to control sample.

Wave length is 500-550 nm. Coloration is stable for 30-40 minutes.

When using the cuvette with different than suggested volume, it is necessary to change proportionally the volumes of components mixture in the ratio:

Serum (distilled water): Reagent A: Reagent B: 0,4 M NaOH = 1:5:5:50

### CALCULATIONS

Calculate the AST activity by calibration curve.

### CONSTRUCTION OF CALIBRATION CURVE

1. For construction of calibration curve prepare samples according to the table:

Solution	Samples numbers				
	1	2	3	4	5
Distilled water, ml	0,1	0,1	0,1	0,1	0,1
Reagent S, ml	-	0,05	0,10	0,15	0,20
Reagent A, ml	0,50	0,45	0,40	0,35	0,30
Reagent B, ml	0,5	0,5	0,5	0,5	0,5
AST activity, mmol/h x l	0	0,5	1,0	1,5	2,0

2. Mix and incubate for 20 minutes at 18-25°C.

3. Add in labeled tubes 5,0 ml of 0,4 M NaOH.

4. Mix and incubate for 10 minutes at 18-25°C.

5. Measure the absorbance of samples 2-5 relative to first sample. Wave-length is 500-550 nm.

6. The AST activity in samples corresponds to table.

7. Construct the curve of absorbance depending on enzyme activity.

1 mmol/h x l = 278 nmol/s x l = 16,7 IU/l

### METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0,36 mmol/h x l = 100 nmol/s x l.

Linearity limit: 4,2 mmol/h x l = 1170 nmol/s x l. For higher values, dilute the sample with 3% solution of albumin in saline solution, and repeat measurement. Multiply the result by dilution coefficient.

Reproducibility (run to run): variation factor < 7 %.

Interference:

Increased concentration of substances, containing ketogroups, cause false increase enzyme activity. Therefore, it is essential for diabetics to test blank sample, where pipette serum after adding Reagent C.

### DIAGNOSTIC CHARACTERISTICS

Aminotransferases catalyze the formation of glutamic acid from 2-oxoglutarate by transfer of amino groups.

AST is normally present in many tissues, but the highest concentrations are found in liver, cardiac muscles, also in kidneys and pancreas.

The serum concentration of AST is elevated in hepatitis and other forms of liver diseases associated with necrosis: infectious mononucleosis, cholestasis, cirrhosis, metastatic carcinoma of the liver delirium tremens, and after administration of various drugs, such as opiates, salicylates or ampicilin.

Also the serum concentration of AST is elevated after cardiac infarction, in disease of skeletal muscles (for example, progressive muscle dystrophy), in acute pancreatitis and other diseases<sup>4,5</sup>.

### BIBLIOGRAPHY

1. V.G. Kolb, V.S. Kamyshnikov. Clinical chemistry reference book., M, Medicine, 1986.
2. Edition by V.V. Menshikov. Laboratory method for clinical assay. Handbook, M, Medicine, 1987.
3. A.M. Goriachkovsky. Reference manual for clinical biochemistry. Odessa, OKFA, 1994.
4. Friedman and Young, Effects of disease on clinical laboratory tests, 3rd ed. AACCC Press. 1987.
5. Tietz NW, Clinical guide to laboratory tests, 2nd ed. Saunders co, 1991

sediment.