# DAC-SpectroMed s.r.l.

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
Cod 2025A150 150 ml  
Cod 2025A600 600 ml

### PRINCIPLE

Aspartate aminotransferase (AST or GOT) catalyzes the transfer of the amino group from aspartate to 2-oxoglutarate, forming oxalacetate and glutamate.

The catalytic concentration is determined from the rate of decrease of NADH color, measured at 340 (334-365) nm, by means of the lactate dehydrogenase (LDH) coupled reaction<sup>1,2,3</sup>.

### CONTENTS AND COMPOSITION

Reagent **A** pH 7,5  
Tris 88 mmol/l  
L-aspartate 260 mmol/l  
Malate dehydrogenase > 900 U/l  
Lactate dehydrogenase > 1500 U/l  
Sodium azide 1,0 g/l  
Reagent **B**  
NADH 240 µmol/l  
2-oxoglutarate 12 mmol/l  
Sodium azide 1,0 g/l

*Harmful! Harmful if swallowed.*

### STORAGE AND STABILITY OF REAGENTS

Reagents at 2-8°C are stable until the expiry date shown on the label. After opening reagents are stable for 30 days.

Working Reagent is stable for 4 weeks at 2-8°C.

Indications of deterioration: Absorbance of the Working Reagent lowers than 1,200 at 334 nm (1 cm cuvette).

### SAMPLES

Serum free of Hemolysis.

AST in serum is stable for 7 days at 2-8°C.

### REFERENCE VALUES

AST/GOT < 37 U/l<sup>4</sup>. These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 340 (334-365) nm.

Cuvettes with 1 cm light path. Stopwatch.

Pipettes for 100 µl, 150 µl, 500 µl and 1,0 ml.

### PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

### REAGENT PREPARATION

Working Reagent:

2 ml Reagent A + 1 ml Reagent B. Mix gently.

### PROCEDURE

Assay conditions

Wavelength: 340 (334-365) nm  
Temperature: 37°C  
Cuvette: 1 cm light path  
Read against: distilled water  
Method: kinetic UV (decreasing)

#### Procedure A

1. Bring the Working Reagent and the photometer to reaction temperature (37°C).
2. Pipette into a cuvette 1 cm light path:  
Working Reagent 1,0 ml  
Sample or Standard 100 µl
3. Mix and insert the cuvette into the photometer. Start the stopwatch.
4. After 1 minute record initial absorbance against the distilled water and record absorbance at 1 minute intervals thereafter for 3 minutes.
5. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ( $\Delta A/\text{min}$ ).

#### Procedure B

1. Bring Reagent A, Reagent B and photometer to reaction temperature (37°C).
2. Pipette into a cuvette 1 cm light path:  
Reagent A 1,0 ml  
Sample or Standard 150 µl
3. Mix and pipette into a cuvette:  
Reagent B 500 µl
4. Mix and insert the cuvette into the photometer. Start the stopwatch.
5. After 1 minute, record initial absorbance against the distilled water and record absorbance at 1 minute intervals thereafter for 3 minutes.
6. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ( $\Delta A/\text{min}$ ).

### CALCULATIONS

The AST/GOT concentration in the sample is calculated using the following general formula:

$$\frac{\Delta A/\text{min}_{\text{Sam}}}{\Delta A/\text{min}_{\text{St}}} \times C_{\text{St}} = C_{\text{Sam}}$$

Calculation using factor

$$340 \text{ nm: Activity (U/l)} = \Delta A/\text{min}_{\text{Sam}} \times 2000$$

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,001  $\Delta A/\text{min}$ =2,00 U/l.

Linearity limit: 0,130  $\Delta A/\text{min}$ =260 U/l.

Repeatability (within run):

Mean Concentration	CV	n
19,1 U/l	1,02 %	20
128 U/l	1,10 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
38,3 U/l	2,07 %	25
134 U/l	1,11 %	25

CV - coefficient of variation n - number of determinations

Interferences: Hemoglobin 1.6 µmol/l (10 g/l), bilirubin 855 µmol/l (0,5 g/l), lipid 3 g/l, glucose 55,5 mmol/l (10 g/l) and ascorbic acid 2,84 mmol/l (0,5 g/l) do not interfere. Other drugs and substances may interfere<sup>4</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

The aminotransferase catalyze the formation of glutamic acid from 2-oxoglutarate by transfer of amino groups.

AST is found in highest concentration in the liver and heart muscle but it is also abundant in skeletal muscle, kidney and pancreas.

The serum concentration of AST is elevated in hepatitis and other forms of hepatic disease associated with necrosis: infectious mononucleosis, cholestasis, cirrhosis, metastatic carcinoma of the liver, delirium tremens, and after administration of various drugs<sup>4,6</sup>.

Serum AST concentration is also elevated after myocardial infarction, in skeletal muscle disease (as progressive muscular dystrophy), in acute pancreatitis or hemolytic disease and other<sup>4,6</sup>.

### BIBLIOGRAPHY

1. Sociedad Española de Química Clínica, Comité científico, Comisión de Enzimas. Método recomendado para la determinación en rutina de la concentración catalítica de la aspartato aminotransferasa en suero sanguíneo humano. Quim. Clin > 1987:6: 235-239.
2. Approved recommendation (1985) on IFCC methods for the Measurement of Catalytic Concentration of Enzymes. Part 2: IFCC method for Aspartate Aminotransferase (EC 2.6.1.1) J Clin Chem Clin Biochem 1986. 24: 497-510.
3. Gelia FJ, Olivella T, Cruz Pastor M, ArenosJ Moreno R, Durban R and Gomez JA. A simple procedure for routine determination of aspartate aminotransferase and alanine aminotransferase with pyridoxal phosphate Clin Chim Acta 1985. 153:241-247.
4. Expert Panel on enzyme of the IFCC, Clin. Chem. Acta, 1976. 70:F19.
5. Young DS, Effects of drugs on clinical laboratory tests, 4th ed. AACCPress, 1995.
6. Friedman and Young, Effects of disease on clinical laboratory tests, 3th ed. AACCPress.1997.