DAC-SpectroMed s.r.l.

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Bili DT – DAC.Lq

TOTAL AND DIRECT BILIRUBIN DSA METHOD with DMSO

For « in vitro » use only
 Store at 2-8°C

Cod 3005B200 100+100 ml
 Cod 3005B1000 500+500 ml

PRINCIPLE

Direct bilirubin in the sample reacts with diazotized sulfanilic acid (DSA) forming a forming a red-colored azocomplex. Dimethylsulfoxide (DMSO) is used as a solvent for determination of total (direct + indirect) bilirubin in the sample. The intensity of coloration, measured at 540(±10) nm, is directly proportional to bilirubin^{1,2} concentration.

CONTENTS AND COMPOSITION

	Code	3005B200	3005B1000
Reagent AD		100 ml	2x250 ml
Sulfanilic acid	30 mmol/l		
Hydrochloric acid	0.2 mol/l		
Reagent AT		100 ml	2x250 ml
Sulfanilic acid	30 mmol/l		
Hydrochloric acid	50 mmol/l		
DMSO	7 mol/l		
Reagent B		8 ml	2x20 ml
Sodium nitrite	29 mmol/l		
Bilirubin standard		-	5 ml

Concentration of the standard is indicated on the vial's label

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.
 Signs of reagent deterioration: presence of suspension, turbidity.

SAMPLES

Serum. Store protected from the light. Bilirubin in serum is stable within 3 months at -20°C, 4 days at 2-8°C, 1 day at 15-25°C.

REFERENCE VALUES

Adults³: Total bilirubin: up to 1.1 mg/dl = 18.8 µmol/l.
 Direct bilirubin: up to 0.2 mg/dl = 3.4 µmol/l.

These values are approximate, it is recommended that each laboratory should establish its own reference values.

ADDITIONAL EQUIPMENT AND REAGENTS

Analyzer, spectrophotometer or photometer with filter 540(±10) nm
 Cuvettes with 1 cm light path. Thermostat for 37°C. Timer. Dispensers for 50 µl, 100 µl and 1,5 ml

PRECAUTION

The kit is only for in vitro use.
 All specimens must be considered potentially hazardous and handled as infectious. Safety precautions for usage of caustic and poisonous substances should be followed.

REAGENT PREPARATION

Reagents are provided ready to use.
 Bilirubin standard: open the vial, avoiding lyophilized material loss; dispense exactly 5 ml of distilled water into the vial. Cap the vial and incubate for 30 minutes at room temperature (15-30°C). Carefully swirl the vial till the complete dissolution on the contents, avoiding the foam forming. Shaking is not permissible!
 Bilirubin in the diluted Bilirubin Standard is stable 4 hours at +25°C; 6 hours at +4°C; 2 weeks at -20°C. Repeated freezing is not permissible!

PROCEDURE

Assay conditions

Method: end point
 Wavelength: 540(±10) nm
 Cuvette: 1 cm light path
 Temperature: 37°C
 Read against: Sample Blank

1. Pipette into labeled test tubes:

Direct bilirubin

	Sample Blank	Sample/Standard
Reagent AD	1,5 ml	1,5 ml
Sample/Standard	100 µl	100 µl
Reagent B	-	50 µl

Total bilirubin

	Sample Blank	Sample/Standard
Reagent AT	1,5 ml	1,5 ml
Sample/Standard	100 µl	100 µl
Reagent B	-	50 µl

NB: Volumes of reagent, standard and samples can be proportionally changed according to the working volumes of cuvettes in the applied analyzers.

2. Mix thoroughly and let the tubes stand for exactly 5 min at 37°C.
 3. Read the absorbance (A_s) of the Sample of the Standard (A_{st}) at 540(±10) nm against the Sample Blank or the Standard Blank.

CALCULATIONS

For the calculation of direct and total bilirubin concentration in sample use the following formula:

$$\frac{A_s}{A_{st}} \times C_{st} = C_s$$

SI UNITS: mg/dl bilirubin = µmol/l bilirubin : 17,1

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0.05 mg/dl = 0.85 µmol/l.

Linearity limit: up to 15 mg/dl = 257 µmol/l.

At higher concentration dilute sample with salina 1:3, repeat the measurement, multiply by 4 (dilution factor)

Reproducibility (within run):

Mean concentration	CV*	n*
0,59 mg/dl = 10,1 µmol/l	3%	20
6,74 mg/dl = 115,2 µmol/l	1%	20

Reproducibility (between run):

Mean concentration	CV*	n*
0,81 mg/dl = 13,8 µmol/l	7,1%	25
1,41 mg/dl = 24,1 µmol/l	6,3%	25

* CV - coefficient of variation n – number of determinations

Interference:

Hemoglobin (10 g/l) do not interfere. Lipemia (triglycerides > 5 g/l) interfere. Other drugs and substances may interfere⁴.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Bilirubin is a waste product derived from the heme moiety of the hemoglobin released from senescent or damaged erythrocytes, that are destroyed in the reticuloendothelial cells. After production, bilirubin is transported to the liver in association with albumin. Inside the hepatocytes bilirubin is conjugated with glucuronic acid and it is excreted into bile. A number of inherited and acquired diseases affect production, uptake, metabolism, and excretion of bilirubin, resulting in hyperbilirubinemia^{3,5}. Conjugated hyperbilirubinemia is associated to a decreased excretion of bile due to liver diseases (hepatitis or cirrhosis) or to intrahepatic or extrahepatic cholestasis.

bilirubin in serum, with use of surfactants as solubilizing agents. Clin Chem 1974; 20: 447-453.

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4. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.

5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

--	--