**DAC-SpectroMed s.r.l.**  
 MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64  
 Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920  
 Email: office@dacspectromed.com  
 www.dacspectromed.com  
 PT MD 11-15796482-001:2003

**Calcium Ars-DAC.Lq**

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|          | Blanc  | Standard | Proba  |
|----------|--------|----------|--------|
| Standard | -      | 10 µl    | -      |
| Proba    | -      | -        | 10 µl  |
| Reagent  | 1,0 ml | 1,0 ml   | 1,0 ml |

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# Calcium Ars-DAC.Lq

CALCIUM

COLORIMETRIC METHOD WITH ARSENAZO III

For « in vitro » use only

Store at 2 - 8°C

Cod 3010C200 200 ml

Cod 3010C500 500 ml

### PRINCIPLE

Calcium with Arsenazo III yields a blue colored complex. The intensity of the colour formed at 620-650 nm is proportional to the calcium concentration in the sample.

### CONTENTS AND COMPOSITION

Reagent pH 6,5  
Arsenazo III 120 mmol/l  
Imidazol Buffer 100 mmol/l  
Calcium Standard 5 ml

Concentration is given on the label.

NB: Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

### STORAGE AND STABILITY OF REAGENTS

Reagent and Calcium Standard are stable at 2-8°C until the expiry date shown on the label when stored tightly closed.

Indications of deterioration: Presence of particles, turbidity, Blank absorbance at 620-650 nm  $\geq 0,05$ .

### SAMPLES

Serum or plasma: Anticoagulants other than heparin should not be used. Calcium is stable 10 days at 2-8°C.

Urine: Collect 24 hours urine specimens in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% V/V). Record the volume. Is stable 10 days at 2-8°C.

Centrifuging and filtering the sample before testing, dilute a sample 1/2 in distilled water. Mix.

### REFERENCE VALUES

Serum or plasma<sup>2</sup>: 8,6 - 10,3 mg/dl = 2,15 - 2,58 mmol/l.

Urine<sup>2</sup>: 100 - 300 mg/ 24 h = 2,5 - 7,5 mmol/ 24 h.

These ranges are given for orientation only.

It is recommended that each laboratory should assign its own reference ranges.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or colorimeter measuring at 620-650 nm. Dropper for 10  $\mu$ l and 1,0 ml.

It is recommended use disposable material to avoid calcium contamination.

### PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

### REAGENT PREPARATION

Reagents are ready for use.

### PROCEDURE

#### Assay conditions

Method: end point  
Wavelength: 620-650 nm  
Temperature: 37°C/15-25°C  
Blank: reagent

- Bring the reagent to the room temperature (16-25)°C
- Pipette into a cuvette:

|          | Blank  | Standard   | Sample     |
|----------|--------|------------|------------|
| Standard | -      | 10 $\mu$ l | -          |
| Sample   | -      | -          | 10 $\mu$ l |
| Reagent  | 1,0 ml | 1,0 ml     | 1,0 ml     |

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

- Mix and incubate for 2 min at 37°C/15-25°C
- Read the absorbance of the Sample and Standard at 620-650 nm, against the Blanc. The color is stable for a least 1 hour.

### CALCULATION

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times K_d = C_{\text{Sample}}$$

### METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0,37 mg/dl = 0,09 mmol/l.

Linearity limit: 30 mg/dl = 7,4 mmol/l.

For higher values dilute sample 1/2 with physiological solution and repeat measurement.

Repeatability (with run):

| Mean concentration       | CV*   | n* |
|--------------------------|-------|----|
| 9,6 mg/dl = 2,4 mmol/l   | 1,7 % | 20 |
| 13,5 mg/dl = 3,38 mmol/l | 1,2 % | 20 |

Repeatability (run to run):

| Mean concentration       | CV*   | n* |
|--------------------------|-------|----|
| 9,6 mg/dl = 2,4 mmol/l   | 2,2 % | 25 |
| 13,5 mg/dl = 3,38 mmol/l | 2,8 % | 25 |

CV\* - coefficient of variation; n\* - number of determination

Sensitivity: 1 mg/dl = 0,018 A.

Interference: Lipemia (triglycerides up to 10 mg/dl) and bilirubin (up to 20 mg/dl) do not interfere.

Hemolysis (hemoglobin above 2,5 g/l) and other drugs and substances may interfere.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

### DIAGNOSTIC CHARACTERISTICS

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones.

A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption.

Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism.

Clinical diagnosis must be made on integration of clinical and laboratory data.

### BIBLIOGRAPHY

- Michaylova V, Ilkova P. Photometric determination of micro amounts of calcium with Arsenazo III Ahal Chim Acta 1971; 53: 194-198.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.