DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
Email: office@dacspectromed.com
www.dacspectromed.com
PT MD 11-15796482-001:2003

Calcium Cr-DAC.Lq

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Calcium Cr-DAC.Lq

CALCIUM

COLORIMETRIC METHOD WITH O-CRESOLPHTALEIN

For «in vitro» use only

Store at 2 - 8°C

Cod 3012C200 200 ml
Cod 3012C500 500 ml

PRINCIPLE

The measurement of calcium in the sample is based on formation of color complex between calcium and o-cresolphthalein in alkaline medium. The intensity of the colour formed at 570 (550-590) nm is proportional to the calcium concentration in the sample.

CONTENTS AND COMPOSITION

Reagent A
Ethanolamine 500 mmol/l
Reagent B
O-Cresolphthalein 0,62 mmol/l
8-Hidroxyquinolein 0,69mmol/l
Calcium Standard 5 ml

Concentration is given on the label.

Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

The reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed.

Indications of deterioration: Presence of particles, turbidity, Blank absorbance at 570 nm $\geq 0,2$.

SAMPLES

Serum or plasma free of haemolysis. Blood anticoagulants with oxalate or EDTA are not acceptable.

Urine: Collect 24 hours urine specimen in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% V/V). Record the volume.

Dilute a sample 1/2 in distilled water. Mix. Multiple results by 2(dilution factor). Calcium is stable 10 days at 2-8°C.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or colorimeter measuring at 570 (550-590) nm. Dropper for 10 μ l and 1,0 ml.

It is recommended use disposable material to avoid calcium contamination.

PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious.

Precautions established for work with caustic and toxic substances should be observed while using the reagents.

WORKING REAGENTS PREPARATION

Reagents are ready for use.

Preparation working reagent: 1 volume Reagent A mix up 1 volume

Reagent B. To make ready directly before using.

Working reagent is stable 1 day at 15-25°C.

PROCEDURE

Method: end point
Wavelength: 570(550-590) nm
Temperature: 37°C/15-25°C
Blank: reagent

1. Pipette into a cuvette:

	Bank	Standard	Sample
Reagent A	0,5 ml	0,5 ml	0,5 ml
Reagent B	0,5 ml	0,5 ml	0,5 ml
Standard	-	10 μ l	-
Sample	-	-	10 μ l

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix and incubate for 5 min at 37°C/15-25°C

4. Read the absorbance of the Sample and Standard at 570 (550-590) nm, against the Blanc.

The color is stable for a least 40 minutes.

CALCULATION

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times K_d = C_{\text{Sample}}$$

REFERENCE VALUES

Serum or plasma²: Adults: 8,5 - 10,5 mg/dl = 2,1 - 2,6 mmol/l.
Children: 10-12 mg/d = 2,5 - 3 mmol/l
Newborns: 8 - 13 mg/dl = 2-3,25 mmol/l

Urine²: Adults: 50 - 300 mg/24 h = 1,25 - 7,5/24 h.
Children: 80-160 mg/24 h = 2 - 4 mmol/24 h.

These ranges are given for orientation only. It is recommended that each laboratory should assign its own reference ranges.

METROLOGICAL CHARACTERISTICS

Sensitivity limit: 0,10 mg/dl = 0,025 mmol/l.

Linearity limit: 15 mg/dl = 3,75 mmol/l.

For higher values dilute sample 1/2 with physiological solution and repeat measurement.

Repeatability(with run):

Mean concentration	CV*	n*
9,08 mg/dl = 2,27 mmol/l	1,97 %	20
15,7 mg/dl = 3,93 mmol/l	1,53 %	20

Repeatability(run to run):

Mean concentration	CV*	n*
9,03 mg/dl = 2,26 mmol/l	1,99 %	20
13,5 mg/dl = 3,57 mmol/l	1,62%	20

CV* - coefficient of variation; n* - number of determination

Sensitivity: 1 mg/dl = 0,034A.

Interference: Lipemia (triglycerides up to 1,25 mg/dl) do not interfere. Substances with calcium and other drugs and substances may interfere.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones. A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption. Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism.

Clinical diagnosis must be made on integration of clinical and laboratory data.

BIBLIOGRAPHY

1. Michaylova V, Ilkova P. Photometric determination of micro amounts of calcium with Arsenazo III Ahal Chim Acta 1971; 53:194-198.
2. Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
3. Yound DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
4. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.