## DAC-SpectroMed s.r.l.

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## Chlor – DAC.Lq

CHLOR

COLORIMETRIC METHOD WITH THIOCYANATE-Hg

For « in vitro » use only


Store at 2 - 8°C

Cod 3016C100

2x50 ml

### PRINCIPLE

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed at 480(440-500) is proportional to the chloride ion concentration in the sample.

### CONTENTS AND COMPOSITION

Reagent	2x50 ml	
Mercuric thiocyanate (II)		2 mmol/l
Ferric nitrate (III)		40 mmol/l
Mercuric nitrate		0,15 mmol/l
Nitric acid		45 mmol/l
Chloride Standard	5 ml	
Chloride aqueous primary standard	125 mmol/l	

*Calibration with the factor or with the aqueous standard may cause bias. It is recommended to calibrate using a serum based standard.*

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label. Chloride standard once open is stable up to 1 month when stored tightly closed at 2-8°C, protect from light and contaminations prevented during their use. Indications of deterioration: Presence of particulate material, turbidity, absorbance of the blank at 480 nm  $\geq$  0,15 (1 cm cuvette).

### SAMPLES

Serum, plasma free of hemolysis.  
CFS, sweat and other body fluids.  
Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.  
Urine 24-hour: dilute a sample ½ in distilled water. Mix. Multiply results by 2. Ion chloride is stable 1 week at room temperature (15-25)°C, in refrigerator (2-8°C) or frozen (-20 °C) temperatures.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 480(440-500) nm. Dropper at 10 µl and 1,0 ml.  
*It is recommended use disposable material to avoid chloride contamination.*

### PRECAUTION

The kit is only for in vitro use.  
However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.  
Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### REAGENT PREPARATION

Reagents are provided ready to use.

### PROCEDURE

Assay conditions

Method: end point  
Wavelength : 480 (440-500) nm  
Light path: 1 cm  
Temperature : 16-25°C/37°C  
Blank: against reagent

1. Bring the Reagent to room temperature (16-25°C).
2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Reagent, ml	1,0	1,0	1,0
Standard, µl	-	10	-
Sample, µl	-	-	10

*NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.*

2. Mix thoroughly and incubate for 5 minute at temperature 37°C or for room temperature at (15-25°C).
3. Read the absorbance of the Standard (A<sub>st</sub>) and the Sample (A<sub>sa</sub>) at 480 (440-500) nm against the Blank.  
The color is stable for a least 30 minutes.

### CALCULATIONS

The chlor concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times K_d = C_{\text{Sample}}$$

### REFERENCE VALUES

Serum, plasma: 95 - 115 mmol/l.

Urine: 110 - 250 mmol/24 h.

CSF: 120 - 130 mmol/l.

Swat: < 60 mmol/l.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### METROLOGICAL CHARACTERISTICS

Detection limit: 1,13 mmol/l.

Linearity limit: 130 mmol/l.

For higher values dilute sample 1+1 with NaCl and multiply result by 2.

#### Repeatability (within run):

Mean Concentration	CV	n
90,7 mmol/l	0,7 %	20
106 mmol/l	0,69 %	20

#### Reproducibility (run to run):

Mean Concentration	CV	n
91,6 mmol/l	0,76 %	20
108 mmol/l	0,74 %	20

CV - coefficient of variation n – number of determinations

Sensitivity: 1 mmol/l = 0,006 A.

Interferences: Bilirubin up to 120 mg/l, bovine serum albumin up to 150 g/l and triglycerides up to 6 g/ did not significantly alter the assay. Hemolysis and other drugs and substances may interfere<sup>4</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure are used.

### DIAGNOSTIC CHARACTERISTICS

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid base balance. Increase in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis. Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis.

### BIBLIOGRAPHY

1. Miller W.G. Chloride. Kaplan A et al. CLIN Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1059-1062.
2. Ibbott F A. et al. New York Academic Press 1965; 101-111.
3. Schoenfeld R G et al. Clin Chem 1964 (10); 533-539.
4. Levinson S S et al. In Faulkner WR et al editors. (9) AACC 1982: 143-148.
5. Yound DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
6. Yound DS. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
7. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC Press, 1999.
8. Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.