## DAC-SpectroMed s.r.l.

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
Cod 3019C200 200 ml

### PRINCIPLE

Very low density lipoproteins (VLDL) and low density lipoproteins (LDL) in the sample precipitate with phosphotungstate and magnesium ions. The supernatant contains high density lipoproteins (HDL). The HDL cholesterol is the spectrophotometrically measured by means of the coupled reactions described below.

### CONTENTS AND COMPOSITION

Reagent A	2x5 ml	
Phosphotungstate		14 mmol/l
Magnesium chloride		2 mmol/l
Reagent B	2x100 ml	pH 6,9
Pipes		90 mmol/l
Phenol		26 mmol/l
Cholesterol esterase		> 1,0 U/ml
Cholesterol oxidase		> 0,3 U/ml
Peroxidase		> 0,65 U/ml
4-aminoantipyrine		0,4 mmol/l
Cholesterol Standard	5 ml	

Cholesterol Standard Aqueous primary standard.  
Calibration with the provided aqueous standard may cause a matrix related bias. It is recommended to calibrate using a serum based standard

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.  
Indications of deterioration:  
Reagent: Presence of particulate material, turbidity, absorbance of the blank over 0,260 at 505 nm (1 cm cuvette).

### SAMPLES

Serum (plasma).  
Cholesterol HDL in serum or plasma is stable for 7 days at 2-8°C.  
Heparin, EDTA, oxalate and fluoride may be used as anticoagulants.

### REFERENCE VALUES

HDL cholesterol concentrations vary considerably with age and sex. The following cut-off point has been recommended for identifying individuals at high risk of coronary artery disease.

HDL Chol:  
Up to 35 mg/dl = 0,91 mmol/l High risk  
> 60 mg/dl = > 1,56 mmol/l Low risk

LDL Chol:  
Up to 150 mg/dl = 3,88 mmol/l High risk  
> 190 mg/dl = > 4,92 mmol/l Low risk

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Thermostatic water bath at 37°C. Desktop centrifuge. Analyzer, spectrophotometer or photometer able to read at 500-550 nm. Doser at 20 µl, 50 µl, 500 µl and 1,0 ml. Taimer.

### REAGENT PREPARATION

Reagents are provided ready to use.

### PRECAUTION

The kit is only for in vitro use.  
However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### PROCEDURE

Assay conditions  
Method: end point  
Wavelength: 500-550 nm  
Light path: 1 cm  
Temperature: 16-25°C/37°C  
Blank: against reagent

### Precipitation

1. Bring the Reagent to room temperature (16-25°C).

2. Pipette into labeled test tubes:

Sample	0,5 ml
Reagent A	50 µl

3. Mix thoroughly and incubate the tubes for 10 minute at room temperature (16-25°C).

4. Centrifuge at a minimum of 400 r.p.m. for 20 minutes.

5. Carefully collect the supernatant.

### Colorimetry

1. Bring the Reagent B to room temperature (16-25°C).

2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Reagent B	1,0 ml	1,0 ml	1,0 ml
Cholesterol Standard	-	20 µl	-
Sample supernatant	-	-	20 µl

3. Mix thoroughly and incubate the tubes for 10 minute at room temperature (16-25°C) or for 5 minutes at 37°C

4. Read the absorbance (A) of the Cholesterol Standard and the Sample at 505 (500-550) nm against the Blank. The color is stable for at least 30 minutes.

### CALCULATIONS

The HDL cholesterol concentration in the sample is calculated using the following general formula:

$$(A_0 / A_{St}) \times C_{St} = C_0$$

mg/dl x 0,0258 = mmol/l

The LDL cholesterol concentration in the sample is calculated using the following general formula:

$$\text{LDL Chol} = \text{Total Chol} - \text{Triglycerides} : 5 - \text{HDL Chol}$$

### QUALITY CONTROL

It is recommended to use the control serum level I and level II to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,3 mg/dl = 0,078 mmol/l.

Linearity limit: 150 mg/dl = 3,9 mmol/l.

Repeatability (within run):

Mean Concentration	CV	n
30 mg/dl = 0,78 mmol/l	3,3 %	20
55 mg/dl = 1,42 mmol/l	2,0 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
30 mg/dl = 0,78 mmol/l	4,2 %	20
55 mg/dl = 1,42 mmol/l	3,2 %	20

Sensitivity: 24 mA x dl/mg = 927 mA x l/mmol

Interferences: Lipemia (triglycerides 10 g/l) does not interfere. Bilirubin (> 10 mg/dl) and hemoglobin (> 5 g/l) may interfere. Other drugs and substances may interfere<sup>4</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure are used.

### DIAGNOSTIC CHARACTERISTICS

HDL play an important part in the removal of cholesterol from tissues and its transportation to the liver for removal as bile acids. Decreased plasma HDL-cholesterol concentrations are positively correlated with the incidence of atherosclerotic diseases, basis of myocardial infarction and cerebrovascular accidents. There are several disease states or environmental influences associated with reduced levels of HDL: acute or chronic hepatocellular diseases, intravenous hyperalimentation, severe malnutrition, diabetes, chronic anemia, myeloproliferative disorders, Tangier disease, analphalipoproteinemia, acute stress, some drugs and smoking.

### BIBLIOGRAPHY

- Naito HK HDL-cholesterol. Kaplan A et al. Clin Chem The C. V. Mosby Co. St Louis. Toronto. Princeton 1984; 1207-1213 and 437.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 2001.
- Tietz NW. Clinical guide to laboratory tests, 3rd ed. AACC, 1995.
- Burtis A et al Tietz Textbook of Laboratory Tests, 3th ed. AACC, 1999.