## DAC-SpectroMed s.r.l.

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
Cod 3020C80 80 ml

### PRINCIPLE

The cholesterol from low density lipoproteins (LDL), very low density lipoproteins (VLDL) and chylomicrons, is broken down by the cholesterol oxidase in an enzymatic accelerated non-color forming reaction.

The HDL cholesterol is the spectrophotometrically measured 600(±20) nm by means of the coupled reactions described below.

### CONTENTS AND COMPOSITION

Reagent A	60 ml	pH 7,0
Good's buffer		
Cholesterin oxidase	> 1 U/ml	
DSBmT	1 mmol/l	
Reagent B	20 ml	pH 7,0
Good's buffer		
Cholesterol esterase	> 1,5 U/ml	
4-aminoantipyrine	1 mmol/l	
Ascorbate oxidase	>3,0KU/ml	
Detergent	2 %	
Peroxidase	> 1 U/ml	
Calibrator HDL/LDL		

Serum. Concentration is given on the label.

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label. Indications of deterioration: Presence of particulate material, turbidity.

Working Calibrator HDL/LDL is stable for 2 week at 2-8°C or for 3 months at -18 °C when frozen in aliquots.

### SAMPLES

Serum (plasma).

Cholesterol HDL in serum or plasma is stable for 7 days at 2-8°C.

Blood anticoagulants with citrat are not acceptable.

### REFERENCE VALUES

HDL cholesterol concentrations vary considerably with age and sex.

The following cut-off point has been recommended for identifying individuals at high risk of coronary artery disease.

Up to 35 mg/dl = 0,91 mmol/l	High risk
> 60 mg/dl = > 1,56 mmol/l	Low risk

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C able to read at 600-700 nm.

Doser at 3, 100, 300 µl. Taimer.

### PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### REAGENT PREPARATION

Reagent A and Reagent B is provided ready to use.

Working Calibrator HDL/LDL: reconstitute with 1,0 ml of distilled water.

### PROCEDURE

Assay conditions

Method:	end point
Wavelength :	600 - 700 nm
Temperature :	37°C

1. Bring the Reagents and the photometer at temperature 37°C.

2. Pipette into labeled test tubes:

	Blanc	Calibrator	Sample
Reagent A	900 µl	900 µl	900 µl
Calibrator HDL/LDL	-	9 µl	-
Sample	-	-	9 µl

3. Mix and incubate for 5 min at 37°C.

4. Read the absorbance of the Sample and Calibrator (A<sub>1</sub>).

5. Pipette into labeled test tubes:

	Blanc	Calibrator	Sample
Reagent B	300 µl	300 µl	300 µl

6. Mix and incubate for 5 min at 37°C.

7. Read the absorbance of the Sample and Calibrator (A<sub>2</sub>) at 600-700 nm against Blanc.

8. Calculate the difference between absorbencies Sample –

$$\Delta A_{\text{Sample}} = (A_2 - A_1)_S \text{ and Calibrator} - \Delta A_{\text{Calibrator}} = (A_2 - A_1)_C$$

### CALCULATIONS

The HDL cholesterol concentration in the sample is calculated using the following general formula:

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Calibrator}}} \times C_{\text{Calibrator}} = C_{\text{sample}}$$

mg/dl x 0,0259 = mmol/l

### QUALITY CONTROL

It is recommended to use the control serum level I and control serum level II to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Detection limit: 2,5 mg/dl = 0,065 mmol/l.

Linearity limit: 200 mg/dl = 5,18 mmol/l.

Repeatability (within run):

Mean Concentration	CV	n
32,9 mg/dl = 0,85 mmol/l	0,8 %	20
50,6 mg/dl = 1,31 mmol/l	0,5 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
32,8 mg/dl = 0,85 mmol/l	1,3 %	40
50,0mg/dl = 1,30 mmol/l	1,5%	40

Interferences: Hemoglobin (10 g/l), lipemia (triglycerides 18 g/l) bilirubin (60mg/dl) does not interfere. Other drugs and substances may interfere<sup>4</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure are used.

### DIAGNOSTIC CHARACTERISTICS

HDL play an important part in the removal of cholesterol from tissues and its transportation to the liver for removal as bile acids. Decreased plasma HDL-cholesterol concentrations are positively correlated with the incidence of atherosclerotic diseases, basis of myocardial infarction and cerebrovascular accidents. There are several disease states or environmental influences associated with reduced levels of HDL: acute or chronic hepatocellular diseases, intravenous hyperalimentation, severe malnutrition, diabetes, chronic anemia, myeloproliferative disorders, Tangier disease, analphalipopro-teinemia, acute stress, some drugs and smoking.

### BIBLIOGRAPHY

1. Naito HK HDL-cholesterol. Kaplan A et al. Clin Chem The C. V. Mosby Co. St Louis. Toronto. Princeton 1984; 1207-1213 and 437.
2. National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung and Blood Institute; 2001.
3. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
4. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 2001.
5. Tietz NW. Clinical guide to laboratory tests, 3rd ed. AACC, 1995.
6. Burtis A et al Tietz Textbook of Laboratory Tests, 3th ed. AACC, 1999.