DAC-SpectroMed s.r.l.

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
Cod 3021C80 80 ml

PRINCIPLE

A specific detergent solubilizes the cholesterol from high density lipoproteins (HDL), very low density lipoproteins (VLDL) and chylomicrons. The cholesterol esters are broken down by cholesterol esterase and cholesterol oxidase in a non-color forming reaction. The LDL cholesterol is then spectrophotometrically measured at 600 (590-700) nm by means of the coupled reactions described below.

CONTENTS AND COMPOSITION

| | | |
|---|-------|-------------|
| Reagent A | 60 ml | pH 7,0 |
| GOOD buffer | | 50 mmol/l |
| Cholesterol esterase | | 380 U/l |
| Cholesterol oxidase | | 380 U/l |
| Catalase | | 400 U/ml |
| N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (TODS) | | 0,45 mmol/l |
| Reagent B | 20 ml | pH 6,3 |
| GOOD buffer | | 50 mmol/l |
| 4-aminoantipyrine | | 1 mmol/l |
| Peroxidase | | 1000 U/l |
| Calibrator HDL/LDL | | |

Serum. Concentration is given on the label.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration: Presence of particulate material, turbidity.

SAMPLES

Serum (plasma).

Cholesterol LDL in serum or plasma is stable for 5 days at 2-8°C.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer thermostatable at 37°C able to read at 600 (590-700) nm.
Thermostatic water bath at 37°C.
Doser at 87, 200, 600 µl. Taimer.

PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

REAGENT PREPARATION

Reagent A and Reagent B is provided ready to use.

Working Calibrator HDL/LDL: reconstitute with 1,0 ml of distilled water. Stable for 2 week at 2-8°C or for 2 months at -18 °C when frozen in aliquots.

PROCEDURE

Assay conditions

Method: end point
Wavelength : 600 (590-700) nm
Temperature : 37°C
Blank: against reagent
1. Bring the Reagents and the photometer at temperature 37°C.
2. Pipette into labeled test tubes:

| | | | |
|--------------------|--------|------------|--------|
| | Blanc | Calibrator | Sample |
| Reagent A | 900 µl | 900 µl | 900 µl |
| Calibrator HDL/LDL | - | 12 µl | - |
| Sample | - | - | 12 µl |

3. Mix and incubate for 5 min at 37°C.
4. Read the absorbance of the Sample and Calibrator (A₁).
5. Pipette into labeled test tubes:

| | | | |
|-----------|--------|------------|--------|
| | Blanc | Calibrator | Sample |
| Reagent B | 300 µl | 300 µl | 300 µl |

6. Mix and incubate for 5 min at 37°C.
7. Read the absorbance of the Sample and Calibrator (A₂) at 600 (590-700) nm against Blanc.
8. Calculate the difference between absorbencies Sample –
 $\Delta A_{\text{Sample}} = (A_2 - A_1)_S$ and Calibrator – $\Delta A_{\text{Calibrator}} = (A_2 - A_1)_C$

CALCULATIONS

The LDL cholesterol concentration in the sample is calculated using the following general formula:

$$\frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Calibrator}}} \times C_{\text{Calibrator}} = C_{\text{Sample}}$$

Conversion factor: mg/dl x 0,02586 = mmol/l

REFERENCE VALUES

The following uniform cut-off points have been established by the US National Cholesterol Education Program and have also been adopted in many other countries for the evaluation of coronary artery disease risk.

| | |
|------------------------------------|-----------------|
| Up to 100 mg/dl = 2,59 mmol/l | Optimal |
| 100-129 mg/dl = 2,59 - 3,34 mmol/l | Near optimal |
| 130-159 mg/dl = 3,37 - 4,12 mmol/l | Borderline High |
| 160-189 mg/dl = 4,14 - 4,90 mmol/l | High |
| > 190 mg/dl = 4,92 mmol/l | Very high |

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the control serum level I and control serum level II to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

METROLOGICAL CHARACTERISTICS

Detection limit: 3,7 mg/dl = 0,096 mmol/l.

Linearity limit: 1000 mg/dl = 25,9 mmol/l.

Repeatability (within run):

| | | |
|-------------------------|-------|----|
| Mean Concentration | CV | n |
| 146 mg/dl = 3,78 mmol/l | 0,7 % | 20 |
| 210 mg/dl = 5,43 mmol/l | 0,6 % | 20 |

Reproducibility (run to run):

| | | |
|-------------------------|-------|----|
| Mean Concentration | CV | n |
| 143 mg/dl = 3,70 mmol/l | 2,0 % | 40 |
| 207 mg/dl = 5,35 mmol/l | 1,7% | 40 |

Interferences: Hemoglobin (5 g/l), lipemia (triglycerides 12,9 g/l) bilirubin (20 mg/dl) does not interfere. Other drugs and substances may interfere⁴.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or a manual procedure are used.

DIAGNOSTIC CHARACTERISTICS

LDL is the main lipoprotein transporting cholesterol from liver to tissues. Increased plasma LDL-cholesterol concentrations are positively correlated with the incidence of atherosclerotic diseases, basis of myocardial infarction and cerebrovascular accidents. There are several disease states or environmental influences associated with increased levels of LDL-cholesterol: nephrosis, diabetes, obesity, some drugs and smoking.

BIBLIOGRAPHY

- Nauck M, Warnick GR, Rifai N. Methods for measurement of LDL-cholesterol: a critical assessment of direct measurement by homogeneous assays versus calculation. Clin Chem 2002; 48: 236-54
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. clinical guide to laboratory tests, 3rd ed. Saunders Co., 1999.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.