DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
 Тел: /+37322/ 574900, 574922/23; факс: /+ 37322/ 574920
 Email: office@dacspectromed.com
 www.dacspectromed.com
 PT MD 11-38623324-002:2002

Cholin - DAC

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Cholin – DAC

Cholinesterase

Kinetic method with Butyrylthiocholine


For « in vitro » use only

Store at 2 - 8°C

Cod 2030C60 20x3 ml
Cod 2030C150 50x3 ml

PRINCIPLE

Cholinesterase hydrolyzed butyrylthiocholine to butyrate and thiocholine. Thiocholine reacts with 5, 5' - dithiobis-2-nitrobenzoic acid (DTNB) to form mercapto-2-nitrobenzoic acid (5-MNBA) according the following reactions:

CONTENTS AND COMPOSITION

Reagent A pH 7.7
Phosphate buffer 50 mmol/l
Reagent B tablets
After dilution
5,5-DTNB 0.25 mmol/l
Butyrylthiocholine 7 mmol/l

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration:

Reagents: Presence of particulate material, turbidity, absorbance of the blank over 1,200 at 405 nm (1 cm cuvette).

SAMPLES

Serum or heparinized plasma collected by standard procedures.
Cholinesterase in serum is stable for 7 days at 2-8°C.

REFERENCE VALUES

Assay temperature	25°C	30°C	37°C
ChE activity (U/l)	3000-9300	3714-11513	4659-14443

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 25, 30 or 37°C and able to read at 405(±10) nm
Cuvettes with 1 cm light path. Doser at 20 µl and 3.0 ml.
General laboratory equipment

PRECAUTION

For in vitro diagnostics only. Handle all patients' samples as potentially dangerous and treat as infectious.
Precautions established for work with caustic and toxic substances should be observed while using the reagents.

QUALITY CONTROL

It is recommended to use the Sera N (cod. 2055S5) and Sera P (cod. 2057S5) to verify the performance of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

REAGENT PREPARATION

Working Reagent:

Dissolve one tablet of Reagent B in 3 ml of Reagent A.

Cap vial and mix gently to dissolve contents.

Stability of Working Reagent – 2 hours at 2-8°C.

PROCEDURE

Method: kinetic
Wavelength: 405(±10) nm
Temperature: 25/30/37°C
Blank: air or distillate water

1. Bring the Working Reagent and the instrument to reaction temperature.

2. Pipette into a cuvette: 25/30°C 37°C

Working Reagent, ml 1.5 1.5

Sample, µl 10

Sample, diluted 1:2 with NaCl 9 g/l, 10 µl

3. Mix and insert the cuvette into the photometer.

4. Record initial absorbance and measure absorbance at 30, 60 and 90 seconds relative to air or distillate water

5. Calculate the difference between consecutive absorbencies, and the average absorbance difference at 30 second (ΔA/30second).

CALCULATIONS

The ChE concentration in the sample (U/l) is calculated using the following general formula:

$$25\text{-}30^\circ\text{C} \quad (\Delta A/30 \text{ sec}) \quad \times \quad 22710 \quad = \quad \text{U/l}$$

$$37^\circ\text{C} \quad (\Delta A/30 \text{ sec}) \quad \times \quad 45420 \quad = \quad \text{U/l}$$

To correct results to other temperatures multiply by temperature conversion factors shown in the following table:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.24	1.55
30°C	0.81	1.00	1.26
37°C	0.64	0.80	1.00

METROLOGICAL CHARACTERISTICS

Detection limit: 7 U/l

Linearity limit: 9084 U/l

For higher values dilute sample 1/5 with distilled water and repeat measurement. Multiply result by 5.

Repeatability (within run):

Mean Concentration CV n
5992 U/l 1,17% 20
3087 U/l 1,82% 20

Reproducibility (run to run):

Mean Concentration CV n
6277 U/l 0,8% 25
3254 U/l 2,03% 25

Sensitivity: 1 U/l = 0.0002 ΔA/30 sec

Interferences:

Haemolysis and hyperbilirubinemia may affect the results. Some drugs and substances may interfere⁴.

These metrological characteristics have been obtained using an analyzer.

Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Cholinesterase is synthesized by the liver and useful for estimation of its function. Low level of cholinesterase indicates failure of protein synthetical function of liver.

Cholinesterase determination is of great importance for identification of patients with atypical form of enzymes and exposure to pesticides organophosphorus.

Cholinesterase activity reduced during acute infections, pulmonary embolism, myodystrophy and myocardial infarct.

Clinical diagnosis should not be made on the findings of a single test

BIBLIOGRAPHY

1. King M, Colineresterase, Kaplan A et al. Clin Chem The C. V. Mosby Co. St Lous. Toronto. Princeton 1984; 1108-1111.
2. Whitaker M, et al Comperason of a Commercially Available Assay System with Two Reference Methods for the Determination of Plasma Collinesterase Variants. , Clin Chem 1983; (29/10); 1746-1760.
3. Young DS, Effects of drugs on Clinical Lab. Tests 4th ed AACC Press, 1995.
4. Young DS, Effects of disease on Clinical Lab. Tests 4th ed AACC Press, 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry 3rd ed. AACC Press, 1997.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC Press, 1995.

result, but should integrate both clinical and laboratory data.