DAC-SpectroMed s.r.l.

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

150 ml

### PRINCIPLE

Creatine kinase (CK) catalyses the phosphorylation of ADP in presence of creatine phosphate forming ATP and creatine. Catalytic concentration of CK is determined by the content of NADPH produced as a result of coupled reactions:

The intensity of produced colour, measured photometrically 340(±10) nm, is proportional to the CK activity<sup>1,2</sup>.

### CONTENTS AND COMPOSITION

	Code	2035C60	2035C150
Reagent A	pH=6.7	50 ml	125 ml
Imidazole buffer	60 mmol/l		
D-Glucose	27 mmol/l		
N-Acetylcysteine	27 mmol/l		
Magnesium acetate	14 mmol/l		
NADP	2,7 mmol/l		
EDTA	2,0 mmol/l		
Hexokinase	> 5000 U/l		
Reagent B		10 ml	25 ml
Imidazole buffer	160 mmol/l		
ADP	11 mmol/l		
AMP	28 mmol/l		
Diadenosine-5-phosphate	55 µmol/l		
EDTA	2 mmol/l		
G-6PDH	> 14000 U/l		
Creatine phosphate	160 mmol/l		

### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

Signs of reagents deterioration: presence of particulate material, turbidity.

### SAMPLES

Serum, plasma. Heparin and EDTA are recommended as anticoagulants. Less than 10% enzyme activity of CK decreases during 1 day at 2-8°C, during 1 hour at 15-25°C.

### REFERENCE VALUES

	Men <sup>1</sup>	Women <sup>1</sup>
37°C	20 - 200 U/l	20 - 180 U/l

These ranges are given for orientation only.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or thermostatic photometer (37°C) able to read at 340 (±10) nm. Pipette for 40 µl, 250 µl, 1,0 ml.

### PRECAUTION

The kit is only for in vitro use. However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### WORKING REAGENT PREPARATION

Mix 5 parts of Reagent A + 1 part of Reagent B. Working reagent is stable 3 weeks at 2-8°C, 48 hours at 15-25°C.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### PROCEDURE

Assay conditions

Method: kinetic (increasing)

Wavelength: 340 (±10) nm

Temperature: 37°C

Read against: reagent

Procedure A:

1. Bring the Working Reagent and the instrument to reaction temperature (37°C)
2. Pipette into a cuvette:

	Blank	Sample
Working reagent	1,0 ml	1,0 ml
Sample	-	40 µl
Distilled water	40 µl	-

3. Mix thoroughly and place cuvette into photometer, start the timer.
4. After 2 minutes measure initial absorbance at 340 nm and absorbances with 1 minute intervals during 3 minutes.
5. Calculate the difference between consecutive absorbances and mean absorbance difference for 1 minute (ΔA/min).

Procedure B:

1. Bring the Reagent A, Reagent B and the photometer to reaction temperature (37°C).
2. Pipette into cuvette with 1 cm light path:

	Blank	Sample
Reagent A	1,0 ml	1,0 ml
Sample	-	40 µl
Distilled water	40 µl	-
Reagent B	200 µl	200 µl

3. Mix thoroughly and pipette into cuvette:
4. Mix thoroughly and place cuvette into photometer, start the timer.
5. After 2 minutes measure initial absorbance at 340 nm and absorbances with 1 minute intervals during 3 minutes.
6. Calculate the difference between consecutive absorbances and mean absorbance difference for 1 minute (ΔA/min).

*NB: Volumes of reagent, standard and samples can be proportionally changed according to the working volumes of cuvettes in the applied analyzers.*

### CALCULATIONS

Enzyme activity of CK is calculated using the following formula:

$$37^\circ\text{C} \quad 340 \text{ nm} \quad \Delta\text{A}/\text{min} \times 4127 = \text{U/l}$$

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,003 ΔA/min=12 U/l.

Linearity limit: 0,436 ΔA/min=1800 U/l.

Repeatability (within run):

Average concentration, U/l	CV%	n
77 U/l	2,5 %	20
624 U/l	1,0 %	20

Reproducibility (run to run):

Average concentration, U/l	CV%	n
83 U/l	2,8 %	20
616 U/l	0,8 %	20

\* CV - coefficient of variation    n – number of determinations

Interferences: glucose up to 7 g/dl, hemoglobin up to 5 g/l and triglycerides up to 7 mmol/l do not interfere with results<sup>1,2</sup>. Other substances and drugs may interfere<sup>3</sup>.

These metrological characteristics were received using analyzer. Results may vary depending on equipment or procedure used.

### DIAGNOSTIC CHARACTERISTICS

Creatine kinase (CK) is an enzyme which is found primarily in skeletal muscle, cardiac muscle, brain tissue and thyroid gland. That's why level of CK in serum is significantly elevated at skeletal muscles, cardiac central nervous system and thyroid disorders.

*In myocardial infarction, CK levels rise during first 3-6 hours and reach the peak values 24 to 36 hours after onset.*

*CK activity returns to the normal level during 3-4 days*

### BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
2. Gerhardt W et al. Creatinine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979; (25/7):1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed. AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC Press, 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3 rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3<sup>rd</sup> AACC 1995.