DAC-SpectroMed s.r.l.

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CK-MB-DAC



CREATINE KINASE MB (CK-MB)
 IMMUNOINHIBITION KINETIC UV METHOD

For « in vitro » use only
 Store at 2-8°C

Code 2037C30 30 ml
 Code 2037C120 120 ml

PRINCIPLE

Specific antibody inhibits both subunits M CK-MM (ck-3) and one subunit M CK-MB (CK-2), that's why subunit B CK-MB is determined (at absence of CK-BB and CK-1)^{1,2} Catalytical concentration of CK-B, corresponding to the half concentration of CK-MB, is determined by NADH concentration, produced as the result coupled reactions of hexokinase (HK) and glucose-6-phosphatedehydrogenase (G6P-DH):

The intensity of produced colour, measured photometrically 340 nm, is proportional to the half concentration of CK-MB.

CONTENTS AND COMPOSITION

	Code	2037C30	2037C120
Reagent A	pH=6.5	25 ml	4x25 ml
Imidazole/MES buffer	100 mmol/l		
D-Glucose	20 mmol/l		
N-Acetylcysteine	20 mmol/l		
Magnesium acetate	10 mmol/l		
NADP	2 mmol/l		
EDTA	2 mmol/l		
ADP	2 mmol/l		
AMP	5 mmol/l		
Hexokinase (HK)	> 2500 U/l		
Anti-human-CK-M	able to inhibit up to 8000 U/l of CK-M		
Reagent B		5 ml	20 ml
P1, P5-di (adenosine-5')	10 µmol/l		
pentaphosphate			
Glucose-6-phosphate dehydrogenase	> 1500 U/l		
Creatine phosphate	30 mmol/l		

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label. Signs of reagents deterioration: presence of particulate material, turbidity.

SAMPLES

Serum or plasma, free of haemolysis. Heparin or EDTA is recommended as anticoagulant.

CK-MB is stable for 2 days at 2-8°C, 1 month at -20°C protected from light.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 340 (±10) nm. Pipette for 40 µl, 1,0 ml.

PRECAUTION

The kit is only for in vitro use.

Patient serum specimens must be handled as potentially dangerous and treated as infectious.

REAGENT PREPARATION

Mix Reagent A and Reagent B in ratio 5+1. Working reagent is stable for 3 weeks at 2-8°C.

QUALITY CONTROL

It is recommended to use the CK-MB control sera (cod 2060K05) to verify the performance of the measurement procedure.

PROCEDURE

Method: kinetic
 Wavelength: 340 nm
 Temperature: 37°C (±1°C)
 Cuvette: 1 cm light path
 Read against: Reagent

1. Bring the Working reagent and the photometer to reaction temperature (37°C).
2. Pipette into cuvette:

	Blank	Sample
Working reagent	1,0 ml	1,0 ml
Sample	-	40 µl
Distilled water	40 µl	-

NB: Volumes of reagent, standard and samples can be proportionally changed according to the working volumes of cuvettes in the applied analyzers

3. Mix thoroughly and place cuvette into photometer, start the timer.
5. After 2 minutes measure initial absorbance at 340 nm and absorbances with 1 minute intervals during 3 minutes against blank.
6. Calculate the difference between consecutive absorbances and mean absorbance difference for 1 minute (ΔA/min).

CALCULATIONS

The CK-MB activity in the sample is calculated using the following general formula:

$$C_{\text{CK-MB}} = \Delta A \times 8352$$

Relative activity index is calculated as follows:

$$\text{CK-MB conc.} / \text{CK Total conc.} \times 100 = \%$$

REFERENCE VALUES

Miocardial infarction incidence is high at the following conditions:

37°C CK-MB > 24 U/l

Relative index of CK-MB: 6-25% of total CK activity.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

Sensitivity: 0,0013 ΔA/min = 11 U/l

Linearity limit: 0,239 ΔA/min = 2000 U/l. At CK-MB activity more than 2000 U/l sample is to be diluted 1:1 by saline and retest. Multiply by 2.

Repeatability (within run):

Mean Concentration	CV, %	n
35,7 U/l	6,31 %	20
132,7 U/l	1,53 %	20

Reproducibility (run to run):

Mean Concentration	CV, %	n
39 U/l	10,4 %	20
136,6 U/l	6,85 %	20

* CV - coefficient of variation n - number of determinations

Interferences: no interferences was observed with hemoglobin up to 1,25 g/dl, bilirubin up to 35 mg/dl, ascorbic acid up to 62 mg/dl, triglycerides up to 250 mg/dl. Presence of high amount of CK-BB or adenylate kinase, macro and micro mitochondrial CK may affect the results⁵. Other drugs and substances may interfere³.

DIAGNOSTIC CHARACTERISTICS

Percent ration of CK-MB activity to total CK activity is usually not more than 6%, but injuries, excessive physical activity, myopathy and myocarditis, stenocardia, circulatory collapse, myocardial infarction, subarachnoid hemorrhage, hypothermia and hyperthermia, malignant hyperpyrexia, muscular dystrophy, Reye's syndrome, some infections, poisoning, lingering arrhythmias, chronic renal insufficiency cause the increase of relative CK-MB index.

Diagnosing of myocardium damage should be based on the clinical observations, CK-MB activity increase value and activity change in time.^{4,7}

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
2. Gerhardt W et al. Creatinine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979; (25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed. AACCC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACCC Press, 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3 rd ed AACCC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd AACCC 1995.