## DAC-SpectroMed s.r.l.

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Cod 3028C50 50 ml

### PRINCIPLE

Copper forms with 4-(3,5-Dibromo-2-pyridylazo)-N-ethyl-N-sulfo-propyl-aniline a chelate complex. The increase of absorbance of this complex can be measured at 580(±10) nm and is proportional to the concentration of total copper in the sample.

### CONTENTS AND COMPOSITION

Reagent	1x50 ml	pH 5.0
Acetate buffer		0.2 mol/l
4-(3,5-dibromo-2-pyridylazo)-N-ethyl-N-sulfo-propylaniline		0.02 mmol/l
Standard	1x5,0 ml	
Copper	200 µg/dl (31.46 µmol/l)	

Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

### STORAGE AND STABILITY OF REAGENTS

The sealed reagents are stable up to the indicated expiry date if stored at 18° - 22°C.

NOTES: If stored at 2°-8°C precipitation may occur. In this case store the reagent at 18-22°C for about 2 hours and mix until the reagent is clear.

### SAMPLES

Serum, plasma and daily urine collected by standard procedures. Copper in serum or plasma is stable for 8 days at 2-8°C. Heparin may be used as anticoagulants. Urine must be filtrate or centrifugate before assay.

### REFERENCE VALUES

Serum, plasma <sup>3</sup> :			
		µg/dl	µmol/l
< 4 months	Sex	8.9 – 46	1.4 – 7.2
4 – 6 months	independent	25 -108	4 - 17
6 months –13 year		51 -121	8 - 19
14 – 19 years	female	70 -159	11 – 25
	male	64 -114	10 - 18
Adults:	female	76 152	12 - 24
	male	70 -140	11 - 22

Urine<sup>3</sup>: 15-70 µg/dl

µg/dl x 0,157 = µmol/l  
 µmol/l x 6,354 = µg/dl

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 580±10 nm. Thermostat at 37 °C. Pipette for 40 µl, 1,0 ml, 2,5 ml.

### PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### REAGENT PREPARATION

Reagent and Standard are provided ready to use.

For daily urine assay Standard must be diluted 1:20 with physiologic saline. Copper content in diluted standard will be 10 µg/dl.

### PROCEDURE

Assay conditions

Method:	end point
Wavelength :	580 (±10) nm
Light path:	1 cm
Temperature :	37°C
Blank:	against reagent

#### For serum, plasma:

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Copper Standard, µl	-	50	-
Sample, µl	-	-	50
Reagent, ml	1.0	1.0	1.0

2. Mix thoroughly and incubate the tubes for 5 minutes at 37°C.

3. Measure the absorbance (A) of the Standard and the Sample at 580 nm against the Blank. The color is stable for at least 30 minutes.

#### For daily urine:

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Distilled water, ml	1.0	-	-
Diluted Copper Standard, µl	-	1.0	-
Sample, µl	-	-	1.0
Reagent, ml	1.5	1.5	1.5

2. Mix thoroughly and incubate the tubes for 5 minutes at 37°C.

3. Measure the absorbance (A) of the Standard and the Sample at 580 nm against the Blank. The color is stable for at least 30 minutes.

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

### CALCULATIONS

The copper concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

For serum (plasma) samples:

If the Copper Standard provided has been used to calibrate (Note 2):

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 200 = C_{\text{Sample}} (\mu\text{g/dl})$$

For daily urine samples:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times 10 = C_{\text{Sample}} (\mu\text{g/dl})$$

For units recalculation the following coefficients may be used:

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Linearity limit: 500 µg/dl = 78.65 µmol/l.

Repeatability (within run): ≤ 4%

Reproducibility (run to run): ≤ 4.5%

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Copper is an essential trace mineral present in all body tissues. Copper, along with iron, helps in the formation of red blood cells. It also helps in keeping the blood vessels, nerves, immune system, and bones healthy. In large amounts, copper is poisonous. A rare inherited disorder Wilson's disease causes deposits of copper in the liver, brain, and other organs. The increased copper in these tissues leads to hepatitis, kidney problems, brain disorders.

Ceruloplasmin and copper tests can help diagnose Wilson's disease. The body normally eliminates excess copper into the bile but also binds some to an enzyme called ceruloplasmin. With Wilson's disease, the binding and excretion processes do not work properly, resulting in decreased concentrations of ceruloplasmin in the blood but increased concentrations of free copper in the blood, urine, and liver.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

Abe A., Yamashita S., Noma A., Clin. Chem., 552-554-35 (1989)  
 Pasquinelli F., Diagnostica e Tecniche di Laboratorio, (pad.:1099-1102) Rossini Ed. (1984)  
 Ciutl R., Galli A. < Giorn. It. Chim. Clin. 12 (2): 91-111 (1987)