DAC-SpectroMed s.r.l.

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
Cod 3042G200 10x20 ml

### PRINCIPLE

Glucose in the sample originates, by means of the coupled reactions described below, forms a complex that can be measured by spectrophotometry.

### CONTENTS AND COMPOSITION

Reagent A	10x20 ml	pH 7.5
TRIS buffer		4 mmol/l
ATP		2,1 mmol/l
Magnesium salt		0,8 mmol/l
Reagent B	10x20 ml	
NAD		2 mmol/l
Hexokinase (HK)		1000 U/l
Glucose-6-phosphate dehydrogenase (G-6-PDH)		1000 U/l

Standard Glucose  
 Aqueous primary standard.

*NB: Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.*

### STORAGE AND STABILITY OF REAGENTS

Reagents and Glucose Standard are stable until the expiry date shown on the label when stored at 2-8°C tightly closed and if contaminations are prevented during their use.

### SAMPLES

Serum free of hemolysis, plasma and urine.  
 Glucose in serum or plasma is stable for 5 days at 2-8°C.  
 Heparin, EDTA, oxalate and fluoride may be used as anticoagulants.

### REFERENCE VALUES

Serum	60-110 mg/dl	3.33-6.11 mmol/l
Urine	0-20 mg/dl	0-1.1 mmol/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 340±10 nm. Thermostatic water bath at 37°C.  
 Pipette for 10 µl, 1,0 ml.

### PRECAUTION

For in vitro diagnostics only.  
 Handle all patients' samples as potentially dangerous and treat as infectious.  
 Precautions established for work with caustic and toxic substances should be observed while using the reagents.

### WORKING REAGENT PREPARATION AND STABILITY

Content of vial Reagent B dissolve in content of vial Reagent A and mix. Working reagent is stable 30 days at 2-8°C, or 7 days at room temperature.

### PROCEDURE

Assay conditions

Wavelength: 340±10 nm  
 Temperature: 16-25/37 °C  
 Cuvette: 1 cm light path  
 Method: endpoint (increasing)  
 Read against: reagent blank

1. Bring the Reagent to room temperature (15-25°C).
2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Glucose Standard, µl	-	10	-
Sample, µl	-	-	10
Reagent, ml	1.0	1.0	1.0

3. Mix thoroughly and incubate the tubes for 10 minute at room temperature (15-25°C) or for 5 minutes at 37°C.

4. Measure the absorbance of tested Sample and Standard relative to Blank. Wave length is 340(±10) nm.

### CALCULATIONS

The glucose concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.  
 Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Detection limit: 0.23 mg/dl = 0.0126 mmol/l  
 Linearity limit: 800 mg/dl = 44.48 mmol/l.

For higher values dilute sample 1/2 with distilled water and repeat measurement.

#### Repeatability (within run):

Mean Concentration	CV, %	n
88 mg/dl = 4.84 mmol/l	1.2	20
326 mg/dl = 17.93 mmol/l	0.9	20

#### Reproducibility (run to run):

88 mg/dl = 4.84 mmol/l	2.7	25
326 mg/dl = 17.93 mmol/l	1.9	25

\* CV - coefficient of variation n – number of determinations

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

Interferences: Hemoglobin (> 5 g/l), lipemia (triglycerides > 1.25 g/l) and bilirubin (10 mg/dl) may interfere. Other substances and drugs may interfere<sup>1</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Glucose is the major source of energy in the body. Insulin, produced by islet cells in the pancreas, facilitates glucose entry into the tissue cells. A deficiency of insulin or a decrease of its effectiveness increases blood glucose. Elevated serum or plasma glucose concentration is found in diabetes mellitus (insulin dependent, non-insulin dependent) and in other conditions and syndromes<sup>2,3</sup>.

Hypoglycemia can occur in response to fasting, or it may be due to drugs, poisons, inborn errors of metabolism or previous gastrectomy<sup>2,5</sup>.

Glucose concentration in arterial blood higher than in venous blood.

For the empty stomach probe 6-8 hours starvation needs.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. Peterson, J.L., Young, D.S., Anal Biochem., 23, (1968); 301
2. Bondar, R.J.L., Mead, D.C. Clin. Chem., 20, (1974); 586
3. Young, D.S., Pestaner, L.C., Gibberman, V., Clin. Chem., 5, 10, (1975)
4. Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997