DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
Email: office@dacspectromed.com
www.dacspectromed.com
PT MD 11-15796482-001:2003

Hb Chr-DAC.Lq

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Hb Chr-DAC.Lq

HEMOGLOBIN

HEMICROMIC METHOD WITH SODIUM LAURYL SULPHATE

For « in vitro » use only

Store at 2-8°C

Cod 3046H3000 3x1000 ml

ADVANTAGES

- Hb Chr don't use cyanide – highly toxic reagent that can be dangerous for the health of laboratory personal, can be reason of allergic reacts and requires additional precautions

- Hb Chr can be used as express-method because of high reaction speed

- Metrological characteristics of Hb Chr are corresponding to the same of HbCN method.

PRINCIPLE

Hemoglobin reacts with sodium lauryl sulphate to form a colored complex - hemichromic.

The intensity of coloration, measured at 540(±20) nm, is proportional to hemoglobin concentration^{1,2}.

CONTENT AND COMPOSITION

Reagent 3 x 10 ml

Concentrate (x100) of transforming reagent

Standard 2 ml

Hemoglobin. Concentration is given on the label.

STORAGE AND STABILITY OF REAGENTS

Reagents at 2-8°C are stable until the expiry date shown on the label when stored tightly closed, contaminations are prevented during their use and are protected from light.

Hemoglobin Standard is stable at 2-8°C 6 months after opening, but not more than expiry date.

SAMPLES

Venous or capillary blood collected by standard procedures.

Hemoglobin in blood is stable for 6 days at 2-8°C.

QUALITY CONTROL

It is recommended to use commercial whole-blood control materials to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 540±20 nm. Timer.

Dropper for 20 µl and 5 ml.

Glass test tubes, not processed EDTA.

PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Hemoglobin Standard provided ready to use.

Working Reagent:

1. Bring the vial with Reagent till 18-25°C. If precipitate present – heat unopened vial up to 50-60°C until precipitate dissolve.

2. For working reagent preparation fresh distilled water must be used.

3. Dilute the contents of the Concentrated Reagent bottle up to 1000 ml with distilled water. Mix thoroughly. Store the diluted reagent in a brown glass bottle. Stable for 6 months at 18-25°C.

Do not freeze it.

PROCEDURE

Assay conditions

Method: end point

Wavelength: 540±20 nm

Temperature: 16-25°C

Blank: distilled water

1. Pipette into labeled test tubes:

	Standard	Sample
Hemoglobin Standard, µl	10	-
Sample, µl	-	10
Working Reagent, ml	2,5	2,5

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

2. Mix thoroughly and incubate the tubes for 5 minutes at room temperature (16-25°C).

3. Measure the absorbance (A) of the Standard and the Sample at 540(±10) nm against the Blank. The color is stable for at least 5 hours.

CALCULATIONS

The hemoglobin concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{sample}}}{A_{\text{Standard}}} \times C_{\text{St}} = C_{\text{Sample}}$$

REFERENCE VALUES

Men³ 130-160 g/l

Women³ 120-140 g/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

Detection limit: not more than 25 g/l.

Linearity limit: 30 - 180 g/l.

Coefficient of variation: less than 2%

Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

Interferences:

Hipertriglicerinemias, leucocytes (leucocytes more than 25x10⁹/l), presence of HbC or HbS, progressive liver diseases, may falsely elevate the hemoglobin concentration. Bilirubin does not interfere. Other drugs and substances may interfere⁴.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Hemoglobin is a red-pigmented, oxygen carrying protein of the red blood cells of vertebrates. It is a tetramer consisting of two pairs of unlike polypeptide chains, carrying an iron-containing porphyrin derivative called heme.

Hemoglobin values are influenced by physiological variations as age, sex, physical exercise, posture, dehydration and altitude, and by pathological processes. Pathological levels are encountered in anemias and in polycythemias. Impaired production of red cells by the bone marrow, excessive blood loss and impaired delivery of red cells to the peripheral blood are the three main causes of anemia^{3,5}.

Increased hemoglobin values are found in polycythemia vera, erythrocytosis, dehydration, newborns, acquired or congenital cyanosis, chronic heart and lung disease, renal cysts and a number of erythropoietin producing tumors^{3,5}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

- Oshiro I., Takenaka T. et al, Clin Biochem 1982, Apr; 15 (2) 83-88.
- Пособие для врачей «Гемихромный метод определения гемоглобина в крови». Москва, 2002.
- Tietz Textbook of Clinical Chemistry/ 3rd edition, Burtis CA. Ashwood ER. WB Saunders Co., 1999.
- Friedman and Young/ Effects of disease on clinical laboratory tests, 3-th ed. AACC Press, 1997.