## DAC-SpectroMed s.r.l.

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PT MD 11-38623324-002:2002

### LDH-DAC. Lq


LACTATE DEHYDROGENASE (LDH)  
KINETIC DGKC METHOD

For « in vitro » use only  
Store at 2-8°C

Cod 2044L50 50 ml  
Cod 2044L200 200 ml

#### PRINCIPLE

LDH catalyses the transformation of Pyruvate to Lactate in the presence of NADH coenzyme. The transformation of NADH to NAD<sup>+</sup> is accompanied by a decrease in absorbance at 340 nm. The change in absorbance correlates with the LDH activity in the serum.

#### CONTENTS AND COMPOSITION

Reagent A	pH 7.5
Phosphate buffer	56 mmol/l
Pyruvate	1.6 mmol/l
Sodium azide	1 g/l
Reagent B	
NADH	240 µmol/l
Sodium azide	1 g/l

Harmful! Harmful if swallowed.

#### STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

Indications of deterioration:

Presence of particulate material, turbidity, absorbance of the blank lower than 1,100 at 340(334-365) nm (1 cm cuvette).

#### SAMPLES

Serum free of hemolysis.

LDH in serum is stable for 24 hours at 2-8°C.

#### REFERENCE VALUES

Serum lactate dehydrogenase activity: 240-480 U/l

Concentrations in newborns and infant are higher than in adults.

Values in men are slightly higher than in women.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

#### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

#### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 340(334-365) nm.

#### PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

#### REAGENT PREPARATION AND STABILITY

One-reagent procedure:

Mix 1 volume of Reagent A with 1 volume of Reagent B

Stability: at 2-8 °C: 4 weeks

#### PROCEDURE

Assay conditions

Wavelength: 340 (334-365) nm

Temperature: 37 °C

Cuvette: 1 cm light path

Read against: distilled water

Method: kinetic (decreasing)

One-reagent procedure:

1. Bring the Working Reagent and the instrument to reaction temperature.

2. Pipette into a cuvette:

Working reagent, ml 1

Sample, µl 25

3. Mix and after 1 minute incubation, measure the change of optical density per minute (ΔA/min) during 3 minutes.

4. Calculate the difference between consecutive absorbencies, and the average absorbance difference per minute (ΔA/min).

Two-reagent procedure:

1. Pipette into a cuvette:

Reagent A, ml 1

Sample, µl 50

2. Mix and after 1 minute incubation, add:

Reagent B, ml 1

3. Mix and after 1 minute incubation, measure the change of optical density per minute (ΔA/min) during 3 minutes.

4. Calculate the difference between consecutive absorbencies, and the average absorbance difference per minute (ΔA/min).

#### CALCULATIONS

Calculation using calibration

$$\frac{\Delta A / \text{min}_{\text{Sample}}}{\Delta A / \text{min}_{\text{Stand}}} \times C_{\text{Stand}} = C_{\text{Sample}}$$

Calculation using factor

334 nm: Activity (U/l) = ΔA/min x 8252

340 nm: Activity (U/l) = ΔA/min x 8450

#### METROLOGICAL CHARACTERISTICS

Detection limit: 0,001 ΔA/min=8.252 U/l

Linearity limit: 0,145 ΔA/min=1200 U/l

For higher values dilute sample 1/10 with distilled water and repeat measurement.

Repeatability (run to run):

Mean Concentration, U/l	CV, %	n
344.6	0.99	25
523.6	3.27	25

Interferences:

Bilirubin (up to 0.5 g/l), lipids up to 10 g/l, glucose up to 10 g/l, ascorbic acid up to 0.5 g/l does not interfere.

Hemolysis may affect the results. Other substances and drugs may interfere<sup>5</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

#### DIAGNOSTIC CHARACTERISTICS

Lactate dehydrogenase (LDH) is present in every cell, it is a tetramere molecule which is a combination of two different tissue components (M-muscle, H-heart).

There are five different isoenzymes:

LDH-1: LDH-2: LDH-3: LDH-4: LDH-5 = 20: 34: 23: 12: 11.

The serum activity is mainly composed of LDH-1, LDH-2 derived from the myocardium and red blood cells, and LDH-5 derived from the liver.

The activities of isoenzymes are different in cases of certain substrates. The inhibitors and pH sensitivities are different.

The various fractions were determined using chromatography in the past but more recently electrophoresis is the method of choice.

The ratio of isoenzymes indicates certain disease states.

The enzyme activity significantly increases 8-12 hours following a myocardial infarction and declines after 4-5 days.

There is an increase in liver diseases, in certain anaemia and tissue injuries.

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3. Young DS. Effects of drugs on clinical laboratory tests, 3rd ed. AACC Press, 1997.

4. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.

5. Friedman and Young. Effects of disease on clinical laboratory tests, 3rd ed. AACC Press, 1997.

6. Ahh. Biol. Clin. 1982; 40:123.