DAC-SpectroMed s.r.l.
MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
Тел: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
Email: office@dacspectromed.com
www.dacspectromed.com
PT MD 11-15796482-001:2003

Magne - DAC.Lq

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Cod 3060M100 100 ml

PRINCIPLE

Magnesium form a purple colored complex when react with calmagite in alkaline solution. Interference by Ca is prevented by the use of EGTA. The intensity of color formed, measured at 500-550 nm, is proportional to the magnesium concentration in the sample.

CONTENTS AND COMPOSITION

Reagent A	1 x 50 ml	
Amino-methyl-propanol		1 mmol/l
EGTA		0.21 mmol/l
Reagent B	1 x 50 ml	
Calmagite		0.30 mmol/l
Magnesium Standard	5 ml	

Aqueous primary magnesium standard.

Concentration is given on the label.

NB: Calibration with the factor or with the aqueous standard may cause bias. It is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Working reagent is stable 4 days at 2-8°C or 24 hours at 15-25°C.

Indications of deterioration:

Reagent: Presence of particulate material, turbidity, absorbance of the blank over 1,400 at 520(±10) nm (1 cm cuvette).

Standard: Presence of particulate material, turbidity.

SAMPLES

Serum free from hemolysis.

Plasma: Heparin recommended as anticoagulant. EDTA and oxalate should not be used as anticoagulants.

Magnesium in serum or plasma is stable for 7 days at 2-8°C.

Urine: 24-hour urine specimen bring to the pH=1 with hydrochloric acid. In case of turbidity appearance hot 10 minutes at 60°C.

Dilute sample 1/10 with distilled water before testing. Multiply result on 10 (Sample Dilution Factor).

Magnesium in urine is stable for 3 days at 2-8°C.

REFERENCE VALUES

Serum and plasma²: 1.6-2.5 mg/dl = 0.66-1.03 mmol/l.

Urine²: 24-244 mg/24-h = 2-21 mmol/24-h.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 500-550 nm. Dropper for 10 µl and 1,0 ml.

PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious.

REAGENT PREPARATION

Mix equal volumes of Reagent A and Reagent B.

QUALITY CONTROL

It is recommended to use the Sera N-DAC (code 2055S5) and Sera P-DAC (code 2057S5) to verify the performance of the measurement procedure.

PROCEDURE

Assay conditions

Method:	end point
Wave length:	520±20 nm
Temperature:	37/15-25 °C
Blank:	reagent

1. Bring the Reagent to the reaction temperature.

2. Pipette into labeled test tubes: (Notes 1, 2)

	Blank	Standard	Sample
Standard, µl	-	10	-
Sample, µl	-	-	10
Reagent, ml	1.0	1.0	1.0

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

2. Mix thoroughly and incubate for 5 minutes at room temperature or 3 minutes at 37 °C.

3. Measure Absorbance (A) of the Standard and the Sample at 520 nm against the Blank.

The colour is stable for at least 30 minutes.

CALCULATIONS

The magnesium concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times \text{Sample Dilution Factor} = C_{\text{Sample}}$$

Conversion factors: mg/dl x 0.412 = mmol/l

0.5 mmol/l = 1 mEq/l = 1.22 mg/dl = 12.2 mg/l

METROLOGICAL CHARACTERISTICS

Detection limit: 0.02 mg/dl = 0.00082 mmol/l.

Linearity limit: 5 mg/dl = 2.05 mmol/l.

For higher values dilute sample 1:2 with distilled water and repeat measurement. Multiply result on 2 (Sample Dilution Factor).

Repeatability (within run):

Mean Concentration	CV	n
2.39 mg/dl	1,18 %	20
4.01 mg/dl	1,73 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
2.27 mg/dl	2,99 %	20
4.14 mg/dl	3,22 %	20

Sensitivity: 1 mg/dl=0.055A

Interferences: Hemolysis and anticoagulants (except heparin), substances and drugs that has magnesium in their content, interfere the results. These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Magnesium is the second more abundant intracellular cation of the human body after potassium, being essential in great number of enzymatic and metabolic processes.

Is a cofactor of all the enzymatic reactions that involve the ATP and is a cofactor of the membrane that maintains the electrical excitability of muscular and nervous cells.

A low magnesium level is found in malabsorption syndrome, diuretic or amiloglucoside therapy, hyperparathyroidism or diabetic acidosis.

Elevated concentration of magnesium is found in uremia, chronic renal failure, glomerulonephritis, Addison's disease or intensive anti acid therapy.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. This reagent may be used in several automated analysers. Instructions for many of them are available on request.

2. It is recommended use disposable material to avoid magnesium contamination.

3. Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC, cod. 2051M5).

BIBLIOGRAPHY

1. Ferrelli E.C. Magnesium/ Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis, Toronto, Princeton1984:1065-1069.

2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.

3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC Press, 2001.

4. Burtis A et al Tietz Textbook of Clinical Chemistry. 3rd ed AACC Press, 1999.

5. Tietz NW, Clinical Guide to Laboratory Tests, 3rd ed AACC Press, 1995.