DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
 Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
 Email: office@dacspectromed.com
 www.dacspectromed.com
 PT MD 11-15796482-001:2003

PhosphoLipid-DAC

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PT MD 11-15796482-001:2003

PhosphoLipid-DAC



QUANTITATIVE DETERMINATION OF PHOSPHOLIPIDS
CHO-POD

For «in vitro» use only
Store at 2-8°C

Cod 3062P50 50 ml

PRINCIPLE

Phospholipids are hydrolysed by phospholipase D and the liberated choline is subsequently oxidized by choline oxidase (CHO) to betaine with the simultaneous production of hydrogen peroxide. In the presence of peroxidase (POD) the hydrogen peroxide couples oxidatively the 4- Aminophenazone (4-AP) and dichlorophenol to form a quinonemine dye:

The intensity of the color measured at 505(±10) nm is proportional to the phospholipids concentration in the sample.

CONTENTS AND COMPOSITION

Reagent A	50 ml	pH 7,55
TRIS		50 mmol/l
Dichlorophenol		2,1 mmol/l
Reagent B	5x10 ml	
PhospholipaseD		400 U/l
Choline oxidase		2200 U/l
Peroxidase		3600 U/l
4- Aminophenazone		1 mmol/l
Phospholipids Standard	5 ml	
Phospholipids aqueous primary standard	300 mg/dl	

STORAGE AND STABILITY OF REAGENTS

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protect from light and contaminations prevented during their use.

Working Reagent is stable 3 weeks at 2-8°C, at 15-25°C – 7 days. Phospholipids Standard once open is stable to 1 month when stored tightly closed at 2-8°C protect from light and contaminations prevented during their use.

Indications of deterioration: Presence of particulate material, turbidity. Blank absorbance at 505 nm ≥0,16.

SAMPLES

Serum or plasma. Phospholipids are stable 3 days at 2-8°C.

REFERENCE VALUES

Adult: 125-275 mg/dl.

The serum phospholipids concentration in normal healthy individuals is in about the same concentration range as total cholesterol. Any change in cholesterol concentration results in a corresponding change in phospholipids in similar direction.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 505(±10)nm. Pipettes for 10 µl and 1,0 ml. Timer. Thermostat at 37°C.

PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

REAGENT PREPARATION

Working Reagents: dissolve the contents of Reagent B in 10 ml Reagent A. Cap and mix gently to dissolve contents. Phospholipids Standard is provided ready to use.

PROCEDURE

Assay conditions

Method:	end point
Wavelength :	505 (490-550) nm
Temperature :	37°C
Blank:	against reagent

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Working Reagent, ml	1.0	1.0	1.0
Phospholipids Standard, µl	-	10	-
Sample, µl	-	-	10

2. Mix thoroughly and incubate for 5 min. at 37°C.

3. Measure the absorbance (A) of the Standard and the Sample at 505(±10) nm against the Blank.

The color is stable for at least 30 minutes.

CALCULATIONS

The phospholipids concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

When utilization Phospholipids Standard:

$$(A_{\text{sample}} / A_{\text{st}}) \times 300 = C_{\text{sample}} \text{ (mg/dl)}$$
$$\text{mg/dl} \times 0,0129 = \text{mmol/l}$$

QUALITY CONTROL

It is recommended to use the Sera N and Sera P (code 2055S5, 2057S5) to verify the performance of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

METROLOGICAL CHARACTERISTICS

Detection limit: 2,54 mg/l

Linearity limit: 600 mg/dl. If the results obtained were greater than linearity limit, dilute the sample ½ with NaCl 9 g/l.

Repeatability (within run):

Mean Concentration	CV	n
121 mg/dl	1,74%	20
221 mg/dl	0,91 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
126 mg/dl	2,31 %	20
225 mg/dl	2,05 %	20

* CV - coefficient of variation n – number of determinations

Sensitivity: 1 mg/dl = 0,0014 A.

Interferences: no influence of ascorbic acid, glucose, bilirubin, uric acid or hemoglobin was found within the range of physiological concentration. Other substances and drugs may interfere³.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Phospholipids are a complex lipid containing phosphorus. Their function as the principal components of cell membranes makes phospholipids essential for all vital processes. The determination of serum phospholipids is an important clinical test in diagnosis of liver diseases, especially obstructive jaundice.

NOTES

1. Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC, cod. 2051M5).
2. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Naito K.N. Lipids. Kaplan A et al. Clin. Chem. The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 918-919 and 570-572.
2. Takeyama M., et al. A new enzymatic method for determination of serum cholinecontaining phospholipids. Clin. Chem. (1977); Acta 79: 93-98.
3. Young DS, Effects of drugs on Clinical Lab. Tests 4th ed AACCC Press, 1995.
4. Young DS, Effects of disease on Clinical Lab. Tests 4th ed AACCC Press, 1995.
4. Burtis A. et al. Tiets Textbook of Clinical Chemistry, 3rd ed. AACCC 1999.
6. Tietz N.W. et al. Clinical Guide to Laboratory Tests, 3rd ed. AACCC 1999.