DAC-SpectroMed s.r.l.

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PT MD 11-15796482-001:2003

Potassium-DAC.Lq

TURBIDIMETRIC DETERMINATION OF POTASSIUM
WITHOUT DEPROTEINIZATION
TETRAPHENYLBORATE REAGENT

For « in vitro » use only
Store at 15-30°C

Cod 3066P50 50 ml
Cod 3066P250 2x125 ml

PRINCIPLE

Potassium ions react with tetraphenylborate sodium and form stable colloid suspension.

The turbidity which has been formed, measured at 578(±10) nm, is proportional to the potassium concentration in the sample.

CONTENTS AND COMPOSITION

Reagent	
Sodium tetraphenylborate	240 mmol/l
Sodium hydroxide	0.5 mmol/l
Stabilizer	
<i>Corrosive!</i> In case of contact with eyes.	
Potassium Standard	5 ml
Potassium	5.0 mmol/l

Aqueous primary standard.

NB: Calibration with the factor or with the aqueous standard may cause bias. It is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label.

SAMPLES

Serum free from hemolysis, plasma

PRECAUTION

The kit is only for in vitro use. All specimens must be considered potentially hazardous and handled as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REFERENCE VALUES

Serum/Plasma	mmol/l	mg/dl
Children	≤ 1 year 3.2-6.1	12.5-23.8
	> 1 year 3.1-5.1	12.1-19.9
Adults	≤ 60 years 3.3-5.1	12.9-19.9
	> 60 years 3.7-5.4	14.5-21.1

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 600±10 nm. Mixer. Dropper for 20 µl and 1,0 ml. Timer.

It is recommended use disposable material to avoid potassium contamination.

REAGENT PREPARATION

Reagent and Standard ready for used.

PROCEDURE

Method:	end point
Wave length:	578±10 nm
Temperature:	25/37 °C
Blank:	reagent

1. Bring the Reagent to the reaction temperature.
2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Standard, µl	-	20	-
Sample, µl	-	-	20
Reagent, ml	1.0	1.0	1.0

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

2. Mix thoroughly and incubate for 5 minutes at room temperature or 3 minutes at 37 °C.
3. To produce a homogeneous turbidity, mix well again before reading.
4. Measure Absorbance (A) of the Standard and the Sample at 578±10 nm against the Blank.

CALCULATIONS

The potassium concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

Conversion factors:

$$\text{mg/dl} \times 0,256 = \text{mmol/l} \quad \text{mmol/l} \times 3,91 = \text{mg/dl}$$

METROLOGICAL CHARACTERISTICS

Linearity limit: 10 mmol/l

For higher values dilute sample 1/2 with physiological solution and repeat measurement. Multiply result on 2 (Sample Dilution Factor).

Coefficient of variation: less than 2%

Interferences:

Hemolysis, turbidity, bilirubin more than 4.0 mg/dl and nitrogen more than 80 mg/dl caused false heightened results.

Some substances and drugs can interfere.

Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Elevated potassium level may be caused by excess potassium consumption, potassium redistribution in organism: exit from cells in extracellular fluid, hemolysis, hard tissue damage, acidosis, dehydration, acute renal insufficiency, Addison's disease, hypofunction of renin-angiotensine-aldosterone system.

Decreased potassium level may be caused by low potassium incoming, loss of potassium with urine, intestinal secrets, sweat, and potassium redistribution in organism: passage from intracellular fluid in cells.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Hillman, G. and Beyer, G: Z. Klin. Chem. U Klin. Biochem. 5:93 (1967).
2. Hoeflmauyr, J.: Praxis und Helfer 8 (1979).
3. Tietz, N.W.: Fundamentals of Clin. Chem. (1976)., 876.