## DAC-SpectroMed s.r.l.

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### PRINCIPLE

Protein in the sample reacts with copper (II) ion in alkaline medium forming a colored complex.

The intensity of coloration, measured at 545(±10) nm, is proportional to total protein concentration<sup>1</sup>.

### CONTENTS AND COMPOSITION

#### Reagent

Sodium hydroxide 0.20 mmol/l  
Copper (II) sulfate 12 mmol/l  
Sodium Potassium tartrate 32 mmol/l  
Potassium iodide 30 mmol/l  
Sodium chloride 0,15 mol/l

*Corrosive! In case of contact with eyes.*

Protein Standard: 5 ml

Bovine albumin = 70 g/l.

Concentration is given on the label.

*Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.*

### STORAGE AND STABILITY OF REAGENTS

Reagents at 2-8°C are stable until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration:

Presence of particulate material, turbidity.

### SAMPLES

Serum (plasma) collected by standard procedures.

Protein in serum (plasma) is stable for 8 days at 2-8°C.

Anticoagulants other than heparin should not be used.

### REFERENCE VALUES

Blood from the umbilical cord: 48 - 80 g/l

Premature: 36 - 60 g/l

Newborn: 46-70 g/l, 1 week: 44-76 g/l, 7 month - 1 year: 51-73 g/l,

1 - 2 year : 56 - 75 g/l, > 3 year: 60 - 80 g/l

Adults: Ambulatory: 64-83 g/l

Recumbent: 60-78 g/l

>60 year: below 2,0 g/l

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 545±10 nm. Dropper for 20 µl and 1,0 ml. Timer.

### PRECAUTION

For in vitro diagnostics only.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

### REAGENT PREPARATION

Reagent and Protein Standard are provided ready to use.

### PROCEDURE

Assay conditions

Method: end point  
Wavelength : 545 (±10) nm  
Light path: 1 cm  
Temperature : 16-25°C  
Blank: against reagent

1. Pipette into labeled test tubes:

	Blank	Standard	Sample
Distilled water, µl	20	-	-
Protein Standard, µl	-	20	-
Sample, µl	-	-	20
Reagent, ml	1,0	1,0	1,0

*NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.*

2. Mix thoroughly and incubate the tubes for 10 minutes at room temperature (16-25°C).

3. Measure the absorbance (A) of the Standard and the Sample at 545(±10) nm. against the Blank.

The colour is stable for at least 1 hours.

### CALCULATIONS

The total protein concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{ST}} \times K_F = C_{\text{Sample}}$$

### METROLOGICAL CHARACTERISTICS

Detection limit: 4.6 g/l

Linearity limit: 150 g/l

For higher values dilute sample 1/2 with distilled water and repeat measurement.

Repeatability (within run):

Mean Concentration	CV	n
44 g/l	1,1 %	20
57 g/l	0,9 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
44 g/l	1,8 %	25
57 g/l	1,9 %	25

Sensitivity: 5 mA x l/g.

Interferences: Hemoglobin (2.5 g/l) and lipemia interfere. Moderate. Bilirubin (20 mg/dl) does not affect the results. Other substances and drugs may interfere<sup>3</sup>.

These metrological characteristics have been obtained using an analyzer.

Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Most of the plasma proteins are synthesized by the liver. The major exception to this is the immunoglobulines which are produced by plasma cells found in the spleen, lymph nodes and bone marrow.

The two general causes of alterations of serum total protein are a change in the volume of plasma water and a change in the concentration of one or more of the serum proteins.

Hyperproteinemia can be caused by dehydration (inadequate water intake, severe vomiting, diarrhea, Addison's disease, diabetic acidosis) or as a result of an increase in the concentration of specific proteins (immunoglobulines in chronic infections, multiple myeloma)<sup>2,4</sup>.

Hypoproteinemia may be caused by hemodilution (salt retention syndromes, massive intravenous infusions), by an impaired synthesis (severe malnutrition, chronic liver disease, intestinal malabsorptive disease), or by an excessive protein loss due to a chronic kidney disease or severe burns<sup>2,4</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### BIBLIOGRAPHY

1. Gornall AG, Bardawill CS, David MM. Determination of serum proteins by means of the Biuret reaction. *J Biol Chem* 1949; 177: 751-766.

2. Tietz NW. Clinical guide to laboratory tests, 3rd ed. Saunders Co, 1999.

3. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCC Press, 1995.

4. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCC Press, 1997.