DAC-SpectroMed s.r.l.

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Sodium - DAC.Lq

SODIUM

PHOTOMETRIC TEST WITH URANYL ACETATE REAGENT

For « in vitro » use only
Store at 18-22°C

Cod 3078S50 50 ml
Cod 3078S250 250 ml

PRINCIPLE

Sodium reacts with uranyl acetate and magnesium and precipitates as the triple salt sodium-magnesium uranyl acetate. The excess of uranyl ions remaining in the solution are forming a yellow-brown colour complex together with thioglycolate acid. The difference between the optical density of the Blank (without sodium precipitation) and the optical density of the specimen reacted, measured at a wavelength of 405 nm, is proportional to the sodium concentration^{1,2}.

CONTENTS AND COMPOSITION

Reagent A
Uranyl acetate 19 mmol/l
Magnesium acetate in ethanol 140 mmol/l

Caustic! Do not pipet by mouth.

Reagent B

Ammonium thioglycolate 550 mmol/l
Ammonia 550 mmol/l

Toxic! Avoid contact with skin and mucous!

Sodium Standard 5 ml

Sodium chloride, sodium concentration is given on the label.

NB Calibration with the provided aqueous standard may cause a matrix related bias. In these cases, it is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable up at 18-22°C until the expiry date shown on the label, stored in the dark in a tightly closed form. At precipitation in Reagent A, it is recommended to warm up with 37°C to dissolve the precipitate.

SPECIMENS

Serum (Plasma)

It is recommended to use anticoagulants that do not contain sodium (lithium, calcium or magnesium heparin).

Sodium is stable for 24 hours at 15-30°C, 2 weeks - at 2-8°C.

REFERENCE VALUES

Serum, plasma: 135-155 mmol/l¹²

These values are for orientation only.

ADDITIONAL EQUIPMENT

Analyzer or photometer with filter able to read at 365/405 nm.

Centrifuge with 4000 rpm. Timer.

Dispensers with a capacity to 20 µl and 1000 µl.

PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious.

PREPARATION REAGENT FOR WORK

The reagents are ready for use.

QUALITY CONTROL

It is recommended to use the serum control pathologic or normal to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme.

PROCEDURE

Method: end point
Wave length: 365/405 nm
Temperature: 25/37 °C
Blank: distilled water

1. Pipette into labeled conical centrifuge tubes with capacity not less than 3 ml:

	Blank	Standard	Test
Reagent A	1,0 ml	1,0 ml	1,0 ml
Bidistilled water	20 µl	-	-
Sodium Standard	-	20 µl	-
Sample	-	-	20 µl

2. Cap and shake all tubes vigorously, incubate for 5 minutes at 16-25°C and again shake vigorously during 60 sec and incubate during 30 min in the dark.

Insufficient mixing or centrifugation lead to false low test results.

3. Centrifuge all tubes at high speed (4000G) for 10 minutes. When removing tubes from centrifuge, take care not to disturb the precipitate.

4. Put in clean labeled tubes supernatant from the tubes with Reagent A and Reagent B:

	Blank	Standard	Test
Supernatant	20 µl	20 µl	20 µl
Reagent B	1,0ml	1,0 ml	1,0 ml

5. After 3 min incubation at 37°C record absorptions of Blank (A_B), Standard (A_{St}) and Test (A_T) at a wavelength of 365/405 nm against distilled water.

CALCULATION

The concentration of sodium (C_T) in the sample is calculated using the following general formula:

$$\frac{A_B - A_T}{A_B - A_{St}} \times C_{St} = C_T$$

Example of Calculation:

Assume the Standard with a sodium value of 150 mmol/l has an absorbance of 0.30 while the Test Specimen and the Blank yielded absorbencies of 0.32 and 0.65 respectively. The sodium concentration of the Test Specimen may then be calculated as follows: (0.65-0.32)/(0.65-0.30) x 150= 141 mmol/l

METROLOGICAL CHARACTERISTICS

Linearity limit: 300 mmol/l

Coefficient of variation: not more than 2,5%

Interferences: Blood calcium, chloride and potassium levels of up to 3 times normal reportedly exert no adverse influence on the procedure. Phosphorus levels exceeding 5 times normal likewise present no problems.

DIAGNOSTIC CHARACTERISTICS

Sodium is the major component of the cations of the extracellular fluid. It is primarily associated with chloride and bicarbonate in the regulation of acid-basic equilibrium.

Sodium's other important function is the maintenance of osmotic pressure of body fluids, thus protecting the body against excessive fluid loss. Sodium also assists in the preservation of normal irritability of muscle and the permeability of the cells.

The main source of body sodium is sodium chloride contained in ingested food. Only about 1/3 of the total body sodium is contained in the skeleton since most of it is contained in the extracellular body fluids.

Sodium is reduced by 1,5-3,0 mmol/l with each rise of glucose in blood for 100 mg/dl.

The combination of hyponatremia with a urine osmolality greater than plasma osmolality, suggests a possible inadequate secretion of ADH.

Na <120 mmol/l can cause weakness, <110 mmol/l - bulbar or pseudobulbar palsy; 90-105 mmol/l cause severe neurological signs and symptoms. Concentrations >155 mmol/l can cause cardiovascular and renal symptoms, especially at lower volume of plasma. Values >160 mmol/l are potentially dangerous for health.

BIBLIOGRAPHY

Guder, W.; Hoffmann G.; Oppitz K. H.: Normalbereiche klin. Chem. Befunde in den Krakenhäusern Munchens (1982)
Trinder, P.: Analzst 76 (1951), 596.