## DAC-SpectroMed s.r.l.

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## TG - DAC.Lq

TRI GLYCERIDES

ENZIMATIC COLORIMETRIC METHOD GPO/PAP

For « in vitro » use only  
Store at 2-8°C

Cod 3085T50 50 ml  
Cod 3085T250 2x125 ml  
Cod 3085T1000 4x250 ml

### PRINCIPLE

Triglycerides in the sample originates, by means of the coupled reactions described below, a colored complex that can be measured by spectrophotometry<sup>1,2</sup>.

Triglycerides + H<sub>2</sub>O  $\xrightarrow{\text{Lipase}}$  Glycerol + Fatty acids

Glycerol + ATP  $\xrightarrow{\text{Glycerolkinase}}$  Glycerol-3-P + ADP

Glycerol-3-P + O<sub>2</sub>  $\xrightarrow{\text{G-3-P-oxidase}}$  Dihydroxyacetone-P + H<sub>2</sub>O<sub>2</sub>

2H<sub>2</sub>O<sub>2</sub> + 4-Aminoantipyrine + 4-Chlorophenol  $\xrightarrow{\text{Peroxidase}}$  Quinoneimine + 4H<sub>2</sub>O

### CONTENTS AND COMPOSITION

**Reagent**  
Pipes 50 mmol/l  
4-chlorophenol 4 mmol/l  
Magnesium chloride 15 mmol/l  
Lipoprotein Lipase > 2 U/ml  
Glycerol kinase > 0,4 U/ml  
Glycerol-3-phosphate oxidase > 1,5 U/ml  
Peroxidase > 2 U/ml  
4-aminoantipyrine 0,4 mmol/l  
ATP 2 mmol/l  
Sodium azide 1 g/l  
Triglycerides Standard 5 ml  
Glycerol equivalent to triolein. Concentration is given on the label. Aqueous primary standard.

**pH 7,2**

### STORAGE AND STABILITY OF REAGENTS

Store at 2-8°C. Reagent and Triglycerides Standard are stable until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Indications of deterioration:

Reagent: Presence of particulate material, turbidity, absorbance of the blank over 0,350 at 492 nm (1 cm cuvette).

Standard: Presence of particulate material, turbidity.

### SAMPLES

Serum collected by standard procedures.

Triglycerides in serum or plasma are stable for 5 days at 2-8°C.

Heparin, EDTA, oxalate and fluoride recommended as anticoagulant.

### REFERENCE VALUES

Men<sup>6</sup>: 0,65-1,85 mmol/l (0,50-1,65 g/l)

Women<sup>6</sup>: 0,55-1,60 mmol/l (0,40-1,40 g/l)

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

### ADDITIONAL EQUIPMENT

Thermostatic water bath at 37°C.

Analyzer, spectrophotometer or photometer able to read at 505 (490-550) nm. Pipettes for 10 µl and 1,0 ml.

### PRECAUTION

The kit is only for in vitro use. However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

### REAGENT PREPARATION

Reagent and Standard are provided ready to use.

### PROCEDURE

Method: end point  
Wavelength : 505 (490-550) nm  
Light path: 1 cm  
Temperature : 37°C  
Blank: against reagent

- Bring the Reagent and the instrument to reaction temperature (37°C).
- Pipette into labeled test tubes:

	Blank	Standard	Sample
Reagent, ml	1,0	1,0	1,0
Distilled water, µl	10	-	-
Triglycerides Standard, µl	-	10	-
Sample, µl	-	-	10

- Mix thoroughly and incubate the tubes for 5 minutes at 37°C.
- Read the absorbance (A) of the Standard and the Sample at 505 (490-550) nm against the Blank.

### CALCULATIONS

The triglycerides concentration (**K<sub>Sample</sub>**) in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{St}}} \times K_{\text{St}} = K_{\text{Sample}}$$

### QUALITY CONTROL

It is recommended to use the Sera N-DAC (code 2055S5) and Sera P-DAC (code 2057S5) to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

### METROLOGICAL CHARACTERISTICS

Detection limit: 0,010 mmol/l.

Linearity limit: 11,4 mmol/l.

Repeatability (within run):

Mean concentration	CV	n
0,8 mmol/l	0,2%	20
10,3 mmol/l	0,65%	20

Reproducibility (run to run):

Mean concentration	CV	n
1,4 mmol/l	2,82%	20
2,4 mmol/l	0,93%	20

Interferences:

Bilirubin (0,5 g/l), glucose (10 g/l) and ascorbic acid (0,1 g/l) do not interfere. Other substances and drugs may interfere<sup>4</sup>.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

### DIAGNOSTIC CHARACTERISTICS

Triglycerides are esters of glycerol and fatty acids coming from the diet or obtained by synthesis mainly in the liver. Triglycerides are transported in plasma by lipoproteins and used by adipose tissue, muscle and other. Their primary function is to provide energy to the cell.

Elevated serum triglycerides levels can be caused by liver disease, diabetes mellitus, nephrosis, hypothyroidism, alcoholism, familial hyperlipoproteinemia IV and V, and other<sup>3,5</sup>.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

### NOTE

Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC, cod. 2051M5).

### BIBLIOGRAPHY

- Bucolo G and David H. Quantitative determination of serum triglycerides by use of enzymes. *Clin Chem* 1973; 19: 476-482.
- Fossati P and Prencipe L. Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. *Clin Chem* 1982; 28: 2077-2080.
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1998.