DAC-SpectroMed s.r.l.

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Cod 3087T500 500 ml
 Cod 3087T1500 1500 ml

PRINCIPLE

Serum β -globulins, γ -globulins and lipoproteins has been precipitated by thymol in TRIS-maleat buffer (pH 7.55). The intensity of turbidity formed, measured at 630(620-660) nm, is proportional to the amount and ratio of individual protein fractions in the sample.

CONTENTS AND COMPOSITION

	3087T500	3087T1500
Reagent T		
при разбавлении pH 7,55	1x11 ml	3x11 ml
Thymol	6,66 mmol/l	
TRIS (oximetyl) aminomethan	10,0 mmol/l	
Maleic acid	3,66 mmol/l	

ATTENTION: Toxic reagent! Do not pipette by mouth!

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 15-25°C until the expiry date shown on the label when stored tightly closed and if contaminations are prevented during their use.

Working reagent is stable 2 months at 15-25°C.

SAMPLES

Serum free from hemolysis.

REFERENCE VALUES

Serum: 0-4 S-H/I

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N and Sera P to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 630 (620-660) nm.

Dropper for 50 μ l and from 3.0 to 5.0 ml.

Thymol-DAC-Standard (cod 3088E125).

PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious.

Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Thymol reagent:

Pour 450 ml of distilled water in 500 ml volumetric flask and put in drop by drop 10 ml of Reagent (dropper tip must be dip into water in volumetric flask). Mix constantly!

Then add distilled water to 500 ml and mix during 10 minutes until dissolved.

PROCEDURE

Assay conditions

Method:	end point
Wave length:	630 (620-660) nm
Temperature:	15-25 °C
Blank:	reagent

1. Pipette into labeled test tubes:

	Blank	Sample
Sample, μ l	-	50
Thymol reagent, ml	3.0	3.0

2. Mix thoroughly and incubate for 30 minutes at 15-25°C.

3. Mix well again before reading.

4. Measure Absorbance (A) of the Sample at 630 nm against the Blank.

5. Calculation must be done by calibration curve using Thymol-DAC-Standard (cod 3088E125).

METROLOGICAL CHARACTERISTICS

Linearity limit: up to 20.0 S-H/I

For higher values dilute sample 1:1 with physiological solution and repeat measurement. Multiply result on 2 (Sample Dilution Factor).

CV: less than 10%

DIAGNOSTIC CHARACTERISTICS

Thymol test intended for liver functional diagnostics.

It is positive in 99-100 % of Hepatitis A, during toxic hepatitis, post-hepatitis and post-necrotic cirrhosis, during collagenous diseases, malaria, viral infections.

Thymol test is negative in 75% cases of obstructive mechanical jaundice and it became positive if complicated with hepatocellular jaundice.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. It is recommended use disposable material to avoid contamination with synthetic washing reagents.

2. Thymol reagent must be used at 20°C.

3. If assayed serum is turbid it is recommended to use 0.05 ml serum in 3 ml of physiological solution as a Blank.

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