DAC-SpectroMed s.r.l.

MD-2012, Moldova, г. Кишинев, ул. Армянская, 47, кв. 64
Тел.: /+37322/ 574900, 574922/23; факс: /+37322/ 574920
Email: office@dacspectromed.com
www.dacspectromed.com
PT MD 11-15796482-001:2003

TIBC - DAC

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PT MD 11-15796482-001:2003

TIBC - DAC.Lq

TOTAL IRON-BINDING CAPACITY (TIBC) OF SERUM
COLORIMETRIC TEST WITH FERROZINE

For « in vitro » use only

Store at 2-8°C

Cod 3083T120 80 ml

PRINCIPLE

Serum proteins are binding ferrous ions in alkaline medium in presence of excess of iron. Then unbinding ferrous ions react with ferrozine forming a colored complex. The intensity of coloration, measured at 560(±10) nm, is proportional to unbinding iron. Difference between known added amount of iron and measured unbinding iron corresponds to iron-binding capacity of serum (IBC). Total iron-binding capacity (TIBC) is a sum of iron-binding capacity and iron content in serum.

CONTENTS AND COMPOSITION

Iron Reagent	20 ml	
Solution of iron salt	89.5 µmol/l	
Buffer Reagent	100 ml	
Tris	0.5 mol/l	
Color Reagent	0.8 ml	
Ferrozine	30 g/l	

STORAGE AND STABILITY OF REAGENTS

Reagents are stable at 2-8°C until the expiry date shown on the label when stored tightly closed.

Indications of deterioration:

Presence of particulate material, turbidity.

SAMPLES

Serum free of hemolysis, plasma

IBC in serum (plasma) is stable for 7 days at 2-8°C.

REFERENCE VALUES

Serum and plasma¹:

Men: 45-75µM

Women: 10-15% less than men

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

QUALITY CONTROL

It is recommended to use the Sera N-DAC (cod. 2055S5) and Sera P-DAC (cod. 2057S5) to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 560±10 nm. Water bath 37°C.

Dropper for 20µl, 500 µl and 2.5 ml.

PRECAUTION

The kit is only for in vitro use.

All specimens must be considered potentially hazardous and handled as infectious.

REAGENT PREPARATION

Reagents are ready for use.

PROCEDURE

Assay conditions

Method: end point

Wave length: 560±10 nm

Temperature: 16-25°C

Blank: distilled water

1. Pipette into labeled test tubes (Notes 1, 2):

	Control	Sample
Buffer Reagent, ml	1.0	1.0
Sample, µl	-	200
Bidistilled Water, µl	200	
Iron Reagent, µl	200	200

2. Mix thoroughly and measure Absorbance (A₁) of sample at 560(±10) nm against the bidistilled water.

3. Add 10 µl Color Reagent, mix thoroughly and let stand the tubes for 10 minutes at 37°C.

4. Measure Absorbance (A₂) of the Sample and Control (A_k) at 560(±10) nm against the bidistilled water.

CALCULATIONS

The iron concentration in the sample is calculated using the following general formula:

$$\frac{A_k + A_1 - A_2}{A_k} \times C_k = C_s$$

METROLOGICAL CHARACTERISTICS

Linearity limit: 0-150 µmol/l

For higher values dilute sample 1:1 with physiological solution and repeat measurement.

Reproducibility: CV < 5%

Interference: Bilirubin does not interfere.

Hemolysis, lipemia, other substances and drugs may interfere.

These metrological characteristics were received using analyzer.

Results may vary depending on equipment or procedure used.

DIAGNOSTIC CHARACTERISTICS

Total iron-binding capacity – it is total amount of iron, binding by serum proteins. The main iron-binding protein is transferrine. Normally only 1/3 of sites in transferrine are occupied, thus serum transferrine has significant reserved iron-binding capacity.

Serum IBC is decreased in hemochromatosis, in acute iron poisoning, in active cirrhosis or acute hepatitis.

Serum IBC is increased in many but not all patients with iron deficiency anemia and in chronic inflammatory disorders, Measurement of serum IBC should not be used as a test for

identification of iron deficiency.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

NOTES

1. Contamination of glassware with calcium will affect the test. Use acid washed glassware or plastic tubes.

2. Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers. In these cases, it is recommended to calibrate using a serum based standard (Multi St-DAC, cod. 2051M5).

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.

2. Young DS. Effects of drugs on clinical laboratory tests, 3rd ed. AACC Press, 1997.

3. Friedman and Young. Effects of disease on clinical laboratory tests, 3rd ed. AACC Press, 1997.