DAC-SpectroMed s.r.l.

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2x50 ml

PRINCIPLE

Unsaturated lipids react with sulphuric acid to form carbonium ions. In a second step the carbonium ions react with phosphovaniline to give a pink colour.

The intensity of the color measured at 520(490-550) nm is proportional to the total lipids concentration in the sample.

CONTENTS AND COMPOSITION

Reagent	2x50 ml	
Phosphovanilline		235 mmol/l
Standard	5 ml	
Total lipids aqueous primary standard		7,5 g/l

NB Calibration with the provided aqueous standard may cause a matrix related bias, especially in some analyzers.

It is recommended to calibrate using a serum based standard

STORAGE AND STABILITY OF REAGENTS

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protect from light and contaminations prevented during their use.

Standard once open is stable to 1 month when stored tightly closed at 2-8°C protect from light and contaminations prevented during their use.

Indications of deterioration:

Presence of particulate material, turbidity. Blank absorbance at 520 nm $\geq 0,32$.

SAMPLES

Serum or plasma. Total lipids are stable 24 h at 15-25°C or 3 days at 2-8°C.

REFERENCE VALUES

Serum or plasma: 450-800 mg/dl.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 520(490-550)nm. Pipettes for 50 μ l, 100 μ l and 1,0, 2,5 ml. Timer. Thermostat at 37°C. Water bath.

Sulfuric acid (H₂SO₄), 98 %, concentrate.

PRECAUTION

The kit is only for in vitro use.

However, all the compounds based on human serum and patient serum specimens must be handled as potentially dangerous and treated as infectious.

REAGENT PREPARATION

Reagents are provided ready to use.

QUALITY CONTROL

It is recommended to use the Sera N (cod. 2055S5) and Sera P (cod. 2057S5) to verify the performance of the measurement procedure. Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do

PROCEDURE

Assay conditions

Method: end point
Wavelength : 520 (490-550) nm
Light path: 1 cm
Temperature : 37°C
Blank: against reagent

1. Workability the sample and Standard with sulfuric acid concentrate and pipette into labeled test tubes:

	Standard	Sample
H ₂ SO ₄ , ml	2.5	2.5
Standard, μ l	100	-
Sample, μ l	-	100

2. Shake thoroughly, incubate for 10 minutes in a boiling water bath at 100°C, after cool in ice water.

3. Transfer into a cuvette.

	Blank	Standard	Sample
Reagent, ml	1.0	1.0	1.0
Standard in H ₂ SO ₄ , μ l	-	50	-
Sample in H ₂ SO ₄ , μ l	-	-	50

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix thoroughly and incubate the tubes for 15 minutes at 37°C.

4. Measure the absorbance (A) of the Standard and the Sample at 520(490-550) nm against the Blank.

The color is stable for at least 1 hour.

CALCULATIONS

The total lipids concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} \times K_d = C_{\text{Sample}}$$

METROLOGICAL CHARACTERISTICS

Detection limit: 0,008 g/l

Linearity limit: 15 g/l. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/l and multiply the result by 2.

Repeatability (within run):

Mean Concentration	CV	n
555mg/dl	2,87 %	20
919 mg/dl	0,7 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
553 mg/dl	1,78 %	25
919 mg/dl	0,63 %	25

* CV - coefficient of variation n - number of determinations

Sensitivity: 1 mg/dl = 0,00066 A.

Interferences: Lipemia (triglycerides < 10g/l) and bilirubin (<20 mg/dl) don't not interfere.

Hemoglobin (>5 g/l) may interfere. Other substances and drugs may interfere³.

These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

The lipids are organic compounds whose more important function is the one to act like fuel. They have an extraordinary yield, favored by possibility of storing itself in remarkable amounts like fatty weave.

Other functions: they are constituent of biological membranes, from protective fatty structures of the internal organs; provide important compounds in the formation with diverse hormones. Great part of the interest in the study of increase of these compounds must to the connection between hyperlipemia and arterosclerosis, diabetes and cardiac disease.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

- Kaplan A et al. Lipids. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 918-919.
- Cottet M.J. et al. Dosage des lipids seriques par la methode sulfo-phosphovanillique (1) de E Chabrol et R. Charonnat. Academie National de Medicine. 1965; 149: 331-338.
- Young DS. Effects of drugs on Clin Lab. Tests. 4th ed AACCC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4 th ed. AACCC Press, 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3 rd ed. AACCC Press, 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3 rd ed. AACCC Press, 1995.

not recover within the acceptable tolerances.