DAC-SpectroMed s.r.l.

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

Cod 3094U100 100 ml
Cod 3094U500 500 ml
Cod 3094U1000 1000 ml

PRINCIPLE

Ammonia and Carbon dioxide (CO₂) are produced when urea is hydrolyzed in presence of Urease. The Ammonia produced in the reaction combines with 2-Oxoglutarate and NADH in the presence of Glutamate dehydrogenase (GLDH) to yield glutamate and NAD⁺.

The NADH/NAD⁺ reaction produces a unique change in absorbance at 340 (334-365) nm, which correlates with the concentration of urea nitrogen in the sample^{1,2}.

CONTENTS AND COMPOSITION

Reagent A	pH 7,6
Tris	100 mmol/l
α-ketoglutarate	9 mmol/l
Urease	> 6500 U/l
GLDH	> 1100 U/l
Sodium azide	9,5 g/l
Reagent B	
NADH	320 μmol/l
Sodium azide	1 g/l
<i>Harmful/Harmful if swallowed.</i>	
Urea Standard	5 ml
urea	10 mmol/l

NB: Calibration with the factor or with the aqueous standard may cause bias. It is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents at 2-8°C are stable until the expiry date shown on the label. Working Reagent is stable for 4 weeks at 2-8°C.

Indications of deterioration:

Absorbance of the Working Reagent lowers than 1,500 at 334 nm (1 cm cuvette).

SAMPLES

Serum free of haemolysis and urine. Urea in serum is stable for 7 days at 2-8°C. Urea in urine is stable for 3 days at room temperature. Before use dilute urine in ratio of 1:100 with distilled water. Do not use anticoagulants containing fluoride or ammonium ions!

REFERENCE VALUES

Serum¹: 2.49-7.47 mmol/l = 0.15-0.40 g/l.

Urine¹: 333-583 mmol/l/24 h = 20-35 g/l/24 h.

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT AND REAGENTS

Analyzer, spectrophotometer or photometer with cell holder thermostatable at 37°C and able to read at 340 (334-365) nm. Pipettes for 10 μl and 1,0 ml. Stopwatch.

PRECAUTIONS

The kit is only for in vitro use.

Handle all patients' samples as potentially dangerous and treat as infectious. Precautions established for work with caustic and toxic substances should be observed while using the reagents.

REAGENT PREPARATION

Working Reagent: 3 ml Reagent A + 1 ml Reagent B. Mix gently.

PROCEDURE

Assay conditions

Method: kinetic
Wavelength : 340 (334-365) nm
Light path: 1 cm
Temperature : 37°C
Blank: distilled water

1. Bring the Working Reagent and the photometer to reaction temperature (37°C).

2. Pipette into a cuvette 1 cm light path:

Working Reagent	1,0 ml
Sample or Standard	10 μl

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix and insert the cuvette into the photometer. Start the stopwatch.

4. After 30 sec record initial absorbance (A₁) against the distilled water and record absorbance for 90 sec (A₂).

CALCULATIONS

The urea concentration in the sample is calculated using the following general formula:

$$\frac{(A_1 - A_2)_{\text{Sam}}}{(A_1 - A_2)_{\text{St}}} \times C_{\text{St}} = C_{\text{Sam}}$$

QUALITY CONTROL

It is recommended to use the Control Serum level I and level II to verify the performance of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

METROLOGICAL CHARACTERISTICS

Detection limit: 0,009 mmol/l = 5.4 mg/dl.

Linearity limit: 66.7 mmol/l = 4 g/l.

Repeatability (within run):

Mean Concentration	CV	n
7,4 mmol/l	1,58 %	20
23,1 mmol/l	1,29 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
7,6 mmol/l	3,26 %	25
25,6 mmol/l	3,51 %	25

* CV - coefficient of variation n - number of determinations

Interferences: Bilirubin 855 μmol/l (0.5 g/l), lipid 10 g/l, glucose 55,5 mmol/l (10 g/l) and ascorbic acid 2,84 mmol/l (0,5 g/l) don't interfere with the assay at the given levels. Other drugs and substances may interfere⁴. These metrological characteristics have been obtained using an analyzer. Results may vary if a different instrument or manual procedures are used.

DIAGNOSTIC CHARACTERISTICS

Urea is synthesized in the liver as a by-product of the deamination of amino acids. Its elimination in the urine represents the major route for nitrogen excretion.

Elevated urea concentration in plasma is found as a result of a high-protein diet, increased protein catabolism, after a gastrointestinal hemorrhage, mild dehydration, shock and heart failure or treatment with glucocorticoids (pre-renal uremia)^{3,5}.

Post-renal uremia is caused by conditions that obstruct urine outflow: nephrolithiasis, tumor or prostatic hypertrophy. The usefulness of urea as an indicator of renal function is limited by the variability of its plasma concentration as a result of nonrenal factors^{3,5}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Talke H, Schubert GE. Enzymatische harnstoffbestimmung in blut und serm im optischen test nach Klinische Wochenschrift 1965; 43:174-175.
2. Gutmann I, Bergmeyer HU. Methods of enzymatic Analysis, ed Bergmeyer HU, Academic Press, NY,1974: 4: 1794-1798.
3. Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
4. Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.