DAC-SpectroMed s.r.l.

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Uric Acid – DAC.Lq

URIC ACID

ENZYMATIC-PHOTOMETRIC URICASE/PAP

For « in vitro » use only

Store at 2-8°C

Code 3096U50 1x50 ml
Code 3096U200 2x100 ml
Code 3096U500 4x125 ml**PRINCIPLE**

Uric acid in the sample transformed by use of uricase and hydrogen peroxide under the catalytic effect of peroxidase (POD) reacts with 4-aminophenazone and 3,5-dichlorophenolsulphonate, forming a red color complex. Intensity of the colouring, measured at 510 (500-550) nm, is proportional to the uric acid concentration.

CONTENTS AND COMPOSITION

Reagent A	pH=8,0
Good's buffer	0.1 mol/l
Ascorbate oxidase	≥ 0.2 KU/l
3,5-dichlorophenol-sulphonate	2,5 mmol/l
Glycerolkinase	≥ 2 KU/l
Reagent B	pH=8,0
Good's buffer	0.1 mol/l
4-aminophenazone	0.8 mmol/l
Peroxidase	≥ 3 KU/l
Uricase	≥ 0.6 KU/l
Uric Acid Standard	5 ml

Aqueous primary standard. Calibration with the factor or with the aqueous standard may cause bias. In these cases, it is recommended to calibrate using a serum based standard.

STORAGE AND STABILITY OF REAGENTS

Reagents are stable until the expiry date shown on the label when stored at 2-8°C. upon opening reagents are stable within 90 days being stored closed in the dark room at 2-8°C.

SAMPLES

Serum, plasma: Uric acid in serum or plasma is stable for 3-5 days at 2-8°C, 6 months at -20°C.

Heparin is used as anticoagulant. EDTA, oxalate and fluoride should not be used as anticoagulants.

Urine: diluted 1:10

REFERENCE VALUES

Serum, plasma:

Men: 3.4-7 mg/dl = 202-416 µmol/l

Women: 2.4-5.7 mg/dl = 143-339 µmol/l

Urine: 250-750 mg/24h = 1.49-4.46 mmol/24 h

These ranges are given for orientation only; each laboratory should establish its own reference ranges.

ADDITIONAL EQUIPMENT

Analyzer, spectrophotometer or photometer able to read at 510 (500-550) nm. Thermostat at 25°C or 37°C.
Dispensers for 25 µl and 1,0 ml.

REAGENT PREPARATION

For Working Reagent preparation mix gently Reagent A and Reagent B in ratio 4+1. Avoid foaming!

Working Reagent is stable for 30 days at 2-8°C.

PRECAUTION

The kit is only for in vitro use.

Patient specimens must be handled as potentially dangerous and treated as infectious.

PROCEDURE

Assay conditions

Method:	end point
Wavelength :	510 (500-550) nm
Temperature :	37°C
Blank:	against reagent

1. Bring the Reagent to room temperature (16-25°C).

2. Pipette into labeled test tubes:

	Blank	Standard	Sample
Working Reagent	1,0 ml	1,0 ml	1,0 ml
Uric acid standard	-	25 µl	-
Sample	-	-	25 µl

NB: Volumes of reagent, standard and samples can be proportionally changed according to the cells working volumes of using analyzers.

3. Mix thoroughly and incubate the tubes for 5 minutes at 37°C.

4. Measure the absorbance (A_{Standard}) of the Standard and of the Sample (A_{Sample}) against the Blank.

CALCULATIONS

The uric acid concentration in the sample is calculated using the following general formula:

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

QUALITY CONTROL

It is recommended to use the Sera N-DAC (cod. 2055S5) and Sera P-DAC (cod. 2057S5) to verify the performance of the measurement procedure.

METROLOGICAL CHARACTERISTICS

Detection limit: 0.2 mg/dl = 11.9 µmol/l.

Linearity limit: 23 mg/dl = 1368 µmol/l.

For higher values dilute sample 1/5 with distilled water and repeat measurement.

Repeatability (within run):

Mean Concentration	CV	n
2.91 mg/dl = 173.4 µmol/l	2.98 %	20
7.39 mg/dl = 440.4 µmol/l	1.56 %	20

Reproducibility (run to run):

Mean Concentration	CV	n
3.14 mg/dl = 187 µmol/l	1.6 %	20
7.83 mg/dl = 466.3 µmol/l	1.42 %	20

* CV - coefficient of variation n - number of determinations

Sensitivity: 33.3 mA x dl/mg = 0.56 mA x l/µmol.

Interferences: Hemolysis (hemoglobin up to 50 mg/l), ascorbic acid up to 20 mg/l, bilirubin up to 20 mg/dl and tricyclerides up to 800 mg/dl do not interfere.

DIAGNOSTIC CHARACTERISTICS

In humans, uric acid is the major product of the catabolism of the purine bases which are obtained partly from the diet and partly from *in vivo* synthesis.

Hyperuricemia is commonly associated with gout, renal failure, leukemia, multiple myeloma, polycythemia, lymphoma, other disseminated tumours, toxicosis of pregnancy, psoriasis, type I glycogenosis, Lesch-Nyhan syndrome, Down's syndrome, polycystic renal disease, chronic lead nephropathy.

Hypouricemia is commonly associated with Wilson's syndrome (hepatolenticular degeneration), Fanconi's syndrome, malignant neoplasms (Hodgkin's disease, multiple myeloma, bronchogenic cancer), xanthinuria, syndrome of pathologic antidiuretic hormonefunction, adenosine deaminase deficiency, purine and nucleoside triphosphatase deficiency; low purine diet.

BIBLIOGRAPHY

- Thefeld C. et al.: Dtsch. Med. Wachr. 98, 380-384 (1973).
- Barham D., Trinder P.: Analyst 97, 142-145 (1972).
- Fossati P., Prencipe L., Bertl G.: Clin. Chem. 26/2, 227-231 (1980).
- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).
- Kaplan L.A., Pesce A. J., ed. Chemistry Theory, Analysis, and Correlation, 3rd ed. St Louis, MO: Mosby, 501-2 (1996).
- Tietz NW. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 624, (1995).